

**United States Department of Agriculture  
Food Safety and Inspection Service, Office of Public Health Science**

CLG-CBX4.00		Page 1 of 12
Title: Screen of Carbadox Metabolite QCA by HPLC-MS-MS		
Revision:NA	Replaces: NA	Effective: 03/28/2011

**Contents**

A.	INTRODUCTION .....	2
B.	EQUIPMENT .....	2
C.	REAGENTS AND SOLUTIONS .....	3
D.	STANDARDS.....	4
E.	SAMPLE PREPARATION.....	5
F.	ANALYTICAL PROCEDURE .....	5
G.	CALCULATIONS .....	7
H.	SAFETY INFORMATION AND PRECAUTIONS .....	8
I.	QUALITY ASSURANCE PLAN .....	9
J.	WORKSHEET.....	11
K.	APPENDIX.....	11
L.	APPROVALS AND AUTHORITIES.....	12

**United States Department of Agriculture  
Food Safety and Inspection Service, Office of Public Health Science**

CLG-CBX4.00	Page 2 of 12	
Title: Screen of Carbadox Metabolite QCA by HPLC-MS-MS		
Revision:NA	Replaces: NA	Effective: 03/28/2011

**A. INTRODUCTION**

Carbadox is an anti-microbial drug that is used as a growth promoter and to prevent swine dysentery and bacterial enteritis.

1. Theory


Homogenized liver tissue is digested with 3M NaOH for 30 minutes in an oil bath. The tissue is acidified, extracted and centrifuged. The extract is then purified by passing through an SCX cartridge. After washing the cartridge with different solvents, 2-quinoxalinecarboxylic acid (QCA), is eluted with 70:30 0.1 M sodium hydroxide: methanol. The eluate is extracted with ethyl acetate, evaporated to dryness, reconstituted with methanol/water, and analyzed by liquid chromatography-tandem mass spectrometry (LC-MS-MS).

Note: 2-Quinoxalinecarboxylic acid is often referred to as quinoxaline-2-carboxylic acid.

2. Applicability

This method is for the screening of QCA, quinoxaline-2-carboxylic acid, the marker for carbadox in swine liver at levels  $\geq 15$  ppb.

3. Structure



Quinoxaline-2-Carboxylic Acid (QCA)

**B. EQUIPMENT**

**Note: Equivalent equipment may be substituted.**

1. Apparatus

- a. Top-Load Balance – sensitive to 0.01 g, Model PM 360, Mettler.
- b. Analytical Balance – sensitive to 0.0001 g, Model AG204, Mettler.
- c. Oil Bath - High Temperature Bath, Fisher Scientific.
- d. Shaker - Eberbach.
- e. Centrifuge - Sorvall RC 4 Centrifuge, Thermo Electron Corporation.

**United States Department of Agriculture  
Food Safety and Inspection Service, Office of Public Health Science**

CLG-CBX4.00	Page 3 of 12	
Title: Screen of Carbadox Metabolite QCA by HPLC-MS-MS		
Revision:NA	Replaces: NA	Effective: 03/28/2011

- f. Nitrogen Evaporator- Model 112 S/N 6831, Organomation.
  - g. Vortex – Genie-2 Cat. No. 12-812, Fisher Scientific.
  - h. Centrifuge tubes – 15 mL glass disposable, with Teflon-lined screw caps, conical bottom, Fisher Scientific.
  - i. Glass centrifuge tubes - 50 mL glass, Heavy Duty, with screw cap, 25 x 150 mm Cat. No. 9826-25, Corning.
  - j. Micropipettes – Adjustable, 10-100 µL, 100-1000 µl, 500-5000 µl, plus tips, Eppendorf.
  - k. Vacuum manifold – SPE 20 port, Waters.
  - l. SPE column – Benzenesulfonic acid cation exchange (SCX) 500-mg/3 mL Cat. No. 2323, Applied Separations.
  - m. Liquid dispensers – Adjustable 1-10 mL, 5-50 mL, Brinkmann.
  - n. Disposable glass culture tubes - 15 mL, 16 x 100 mm Cat. No. 73500 16100, Kimble.
  - o. pH meter – Ultrabasic pH meter, UB-5, Denver Instruments.
2. Instrumentation
- a. Waters HPLC/MS/MS – Waters Quattro Micromass, Alliance Column – Pursuit XRs C18 µm 100 x 2 mm

**C. REAGENTS AND SOLUTIONS**

**Note: Equivalent reagents and solutions may be substituted.**

1. Reagents
- a. Sodium hydroxide - Pellets, Cat. No. SX0590-14, EMD.
  - b. 3N Sodium Hydroxide - Cat. No. BDH 3472-1, BDH.
  - c. 1N Sodium Hydroxide - Cat. No. 7450-32, Ricca Chem Company.
  - d. Hydrochloric acid – Reagent grade, Cat. No. 9535-33, J.T Baker.
  - e. 1N Hydrochloric acid - Cat. No. BDH 3201-2, BDH.
  - f. Ethyl Acetate - Reagent grade, Cat. No. 100-4, Burdick and Jackson.
  - g. Sodium phosphate monobasic - Cat. No. S8282, Sigma Aldrich.
  - h. Sodium phosphate dibasic – Cat. No. S7907, Sigma Aldrich.
  - i. Methanol – High Purity Solvent, Cat. No. 230-4, Burdick and Jackson.

**United States Department of Agriculture  
Food Safety and Inspection Service, Office of Public Health Science**

CLG-CBX4.00	Page 4 of 12	
Title: Screen of Carbadox Metabolite QCA by HPLC-MS-MS		
Revision:NA	Replaces: NA	Effective: 03/28/2011

- j. Nitrogen – Airgas.
- k. Argon – Airgas.
- l. Formic Acid - Cat. No 94318, Fluka Analytical.

2. Solutions

- a. 3 M NaOH:  
Add 120 g of NaOH pellets and 1000 mL of Milli Q water and mix.
- b. 0.2 M sodium phosphate monobasic:  
Weigh out 13.93 g of sodium phosphate monobasic. Then add 500 mL of Milli Q water.
- c. 0.2 M sodium phosphate dibasic:  
Weigh out 53.65 g of sodium phosphate dibasic. Add 500 mL of Milli Q water.
- d. 0.1 M sodium phosphate, pH 8-8.5:  
Add 15.9 mL of 0.2 M sodium phosphate monobasic and 285 mL of 0.2 M sodium phosphate dibasic. Add 600 mL of Milli Q water. Measure pH about 8-8.5.
- e. 0.1 M NaOH:  
Add 100 mL of 1 N NaOH and bring to volume in a 1000 mL class A volumetric flask using Milli Q water.
- f. 0.1 M NaOH:MeOH (70:30):  
Add 700 mL of 0.1 M NaOH to 300 mL of methanol.
- g. 0.1 M HCl:  
Add 100 mL of 1 N HCl and bring to volume in a 1000 mL class A volumetric flask using Milli Q water.
- h. MeOH: H<sub>2</sub>O (10:90):  
Add 10 mL of methanol to 90 mL of Milli Q water and mix.
- i. 0.3 % formic acid:  
Dilute 3 mL of formic acid to 1000 mL with Milli Q water in a 1 L volumetric flask.

**D. STANDARDS**

**Note: Equivalent standards may be substituted.**

**United States Department of Agriculture  
Food Safety and Inspection Service, Office of Public Health Science**

CLG-CBX4.00	Page 5 of 12	
Title: Screen of Carbadox Metabolite QCA by HPLC-MS-MS		
Revision: NA	Replaces: NA	Effective: 03/28/2011

1. Source
  - a. 2-Quinoxalinecarboxylic acid (15 µg/mL in methanol) (QCA) – Cat. No. 91819 Absolute Standards, Aldrich Chemicals.
2. Preparation
  - a. QCA solution (0.15 µg/mL):

Add 250 µL of 15 µg/mL QCA into a 25 mL class A volumetric flask. Dilute to volume using methanol. Mix thoroughly and store at 2 - 4 °C. Prepare every six months.

**E. SAMPLE PREPARATION**

Homogenize liver samples using a blender or food processor.

**F. ANALYTICAL PROCEDURE**

1. Extraction
  - a. Weigh out 5.00 ± 0.10 g into a 50 mL glass tube. Weigh out three blank samples. One will be used as a blank, another as a positive control and the third as a blind check sample.
  - b. Fortify with 500 µL of 0.15 µg/mL QCA (equivalent to 15 ppb recovery) and vortex.
  - c. Let stand for 10 min.
  - d. Add 10 mL of 3 M NaOH to each tube and vortex.
  - e. Place tubes in an oil bath that has been heated to 95- 100 °C. Leave tubes in the oil bath for 30 minutes. After 30 minutes remove the tubes from the oil bath and cool them to room temperature.

Stopping point: The samples may be stored in the refrigerator for 24 hours.
  - f. Once the tubes are cooled to room temperature add 4 mL of concentrated HCl. Cap and vortex for 30 seconds.

Caution: Reaction is highly exothermic. Take care to work in a hood and wear eye and hand protection.
  - g. Add 6 mL of ethyl acetate and shake for one minute. Centrifuge at 2000 g at 4 °C for 10 minutes.
  - h. Transfer the upper layer to a 50 mL polypropylene tube.
  - i. Repeat step g. and combine the extracts.
  - j. Add 8 mL of 0.1 M sodium phosphate solution (C.2.d.) and shake for 1 minute on

**United States Department of Agriculture  
Food Safety and Inspection Service, Office of Public Health Science**

CLG-CBX4.00	Page 6 of 12	
Title: Screen of Carbadox Metabolite QCA by HPLC-MS-MS		
Revision:NA	Replaces: NA	Effective: 03/28/2011

a shaker. Centrifuge at 2000 g at 4 °C for 10 minutes.

- k. Aspirate the upper organic layer.
- l. Take a 4 mL aliquot and add it to a 10 mL glass tube containing 1 mL of concentrated, HCl.

2. Solid Phase Extraction

- a. Prime the SCX SPE cartridge with 5 mL of methanol, followed by 5 mL 0.1 M HCl.
- b. Add 2 mL of the sample.
- c. Wash the cartridge with 5 mL of 0.1 M HCl.
- d. Elute with 3 mL of 0.1M NaOH:MeOH (70:30) into a 15 mL centrifuge tube.
- e. Add 300 µL of concentrated HCl to each sample.
- f. Extract with 2 mL of ethyl acetate. Vortex for 15 seconds. Centrifuge at 2000 g at 4 °C for 10 minutes. Collect the upper organic layer into a 10 mL test tube.
- g. Repeat step f. two more times and combine the extracts.
- h. Evaporate the combined ethyl acetate extracts at 50-60 °C to dryness.
- i. Add 100 µL of MeOH:H<sub>2</sub>O (10:90) and vortex for 15 seconds.
- j. Transfer to a vial containing an insert.  
Stopping point: Samples may be left for up to 48 hours at room temperature.
- k. Analyze on the HPLC/MS/MS.

3. Instrumental settings

- a. Set HPLC Parameters:  
Column Temperature: 30°C  
Injection Volume: 15 µL  
Total column flow: 0.35 mL/min  
Pump gradient: (Total run time of 11.01 min)  
Mobile Phase A: Acetonitrile  
Mobile Phase B: 0.3 % formic acid  
Mobile Phase C: Methanol

Gradient set of Mobile phases as follows:

**United States Department of Agriculture  
Food Safety and Inspection Service, Office of Public Health Science**

CLG-CBX4.00	Page 7 of 12	
Title: Screen of Carbadox Metabolite QCA by HPLC-MS-MS		
Revision:NA	Replaces: NA	Effective: 03/28/2011

Time	A%	B%	C%
0.00	2	90	8
10.00	10	50	40
10.01	16	20	64
11.01	2	90	8

- b. Set Mass Spectrometric Parameters. System may be adjusted to ensure optimum performance.

Capillary Voltage 3.00  
Cone Voltage 20.00  
Extractor Voltage 3.00  
RF Lens Voltage 0.3  
Source Temperature 125 °C  
Desolvation Temperature 450 °C

- c. Measure the detector response (peak area) for the product ion transitions  
175 → 129, Collision energy 30 eV, Dwell time 0.3  
175 → 102, Collision energy 18 eV, Dwell time 0.3  
175 → 131, Collision energy 14 eV, Dwell time 0.3

4. Chromatograms

See Section K.1.

**G. CALCULATIONS**

1. Screening Analysis for QCA

- a. A sample will be considered positive if:
- The retention time of the analyte in the sample matches the retention time of either the external standard or the positive control run under the same experimental conditions to within  $\pm 5\%$ .
  - All product ion transitions listed in Section F.3.c. are present.

**United States Department of Agriculture  
Food Safety and Inspection Service, Office of Public Health Science**

CLG-CBX4.00	Page 8 of 12
Title: Screen of Carbadox Metabolite QCA by HPLC-MS-MS	
Revision:NA	Replaces: NA
Effective: 03/28/2011	

- iii. Peaks for all transitions have a signal-to-noise (s/n) ratio  $\geq 3$ .
- iv. The most abundant product ion area is  $> 25\%$  of the positive control.

**H. SAFETY INFORMATION AND PRECAUTIONS**

- 1. Required Protective Equipment — Safety glasses, disposable gloves, lab coat.
- 2. Hazards

<i>Procedure Step</i>	<i>Hazard</i>	<i>Recommended Safe Procedures</i>
Methanol, Ethyl Acetate	Flammable, poisonous, inhalation will cause headaches, fatigue, nausea.	Wear gloves and work in the hood. Use protective eyewear. Avoid contact with skin, eyes.
Hydrochloric acid, Formic acid	Corrosive	Wear PPE, and avoid contact with skin.
QCA	Carbadox has been shown to cause cancer in laboratory animals but when fed to swine, is metabolized or transformed over a relatively short period of time.  May cause photosensitization. This substance is possibly carcinogenic to humans.	Wear PPE, avoid skin contact.  Exercise appropriate precautions to minimize direct contact with skin or eyes and prevent inhalation.
Sodium Hydroxide	Corrosive	Wear PPE, and avoid skin contact.
Sodium Phosphate	Slight reactivity rating. Irritant due to its slight acidic nature. Chronic exposure may result in calcium phosphate deposits in the kidneys.	Wear gloves and work in the hood. Use protective eyewear. Avoid contact with skin and eyes.

- 3. Disposal Procedures

<i>Procedure Step</i>	<i>Hazard</i>	<i>Recommended Safe Procedures</i>
Methanol, Ethyl Acetate	Flammable, Poisonous, inhalation will cause headache, fatigue nausea.	Store waste in a tightly sealed container away from non-compatibilities in a cool, well ventilated, flammable liquid

**United States Department of Agriculture  
Food Safety and Inspection Service, Office of Public Health Science**

CLG-CBX4.00	Page 9 of 12	
Title: Screen of Carbadox Metabolite QCA by HPLC-MS-MS		
Revision:NA	Replaces: NA	Effective: 03/28/2011

		storage area/cabinet for disposal in accordance with local, state, and Federal regulations.
Hydrochloric Acid, Formic acid	Corrosive. The neutralization is an exothermic reaction which may spatter.	Collect waste and store in a tightly sealed container. Store away from non-compatibles in a cool, well ventilated, acid liquid storage area/cabinet for disposal in accordance with local, state and Federal regulations or neutralize and dispose in accordance with local state, and Federal regulations.
QCA	Carbadox has been shown to cause cancer in laboratory animals but when fed to swine, is metabolized or transformed over a relatively short period of time.  May cause photosensitization. This substance is possibly carcinogenic to humans.	Disposal in accordance with local, state, and Federal regulations.

**I. QUALITY ASSURANCE PLAN**

1. Performance Standard

- a. The external standard and the positive control must be positive according to the criteria listed in G.1.
- b. The area response of the negative control must be  $\leq 5\%$  of the positive control.

2. Critical Control Points and Specifications

<i>Record</i>	<i>Acceptable Control</i>
a. Fortify weighed control sample	Let the sample stand for 10 minutes.
b. Sample weight	5.00 ± 0.10 g

3. Readiness To Perform

**United States Department of Agriculture  
Food Safety and Inspection Service, Office of Public Health Science**

CLG-CBX4.00	Page 10 of 12	
Title: Screen of Carbadox Metabolite QCA by HPLC-MS-MS		
Revision:NA	Replaces: NA	Effective: 03/28/2011

- a. Familiarization
    - i. Phase I: Standards- Analyze a 0.15 µg/mL external standard on three different days.
    - ii. Phase II: Fortified samples- Analyze at least 4 fortified tissues per day over 3 different days. They shall be fortified at 15 ppb.  
NOTE: Phase I and Phase II may be performed concurrently.
    - iii. Phase III: Check samples for analyst accreditation.
      - (a) A total of 8 blind check samples, consisting of one or two unfortified blanks and the remainder of the samples fortified at 15 ppb.
      - (b) Report analytical findings to Supervisor and Quality Assurance Manager (QAM).
      - (c) Authorization from QAM and Supervisor is required to commence official analysis.
  - b. Acceptability criteria.  
Refer to I. 1.
4. Intralaboratory Check Samples
- a. System, minimum contents.
    - i. Frequency: One per week per analyst when samples analyzed
    - ii. Records are to be maintained
  - b. Acceptability criteria.  
Refer to I. 1.  
If unacceptable values are obtained, then:
    - i. Stop all official analyses by that analyst for this method.
    - ii. Take corrective action.
5. Sample Acceptability and Stability
- a. Martix: Swine liver
  - b. Condition: Cold
  - c. Sample storage:
    - i. Time: 6 months
    - ii. Condition: Frozen at -20 °C.

**United States Department of Agriculture  
Food Safety and Inspection Service, Office of Public Health Science**

CLG-CBX4.00	Page 11 of 12	
Title: Screen of Carbadox Metabolite QCA by HPLC-MS-MS		
Revision:NA	Replaces: NA	Effective: 03/28/2011




6. Sample Set
  - a. External Standard
  - b. Negative Control (blank tissue)
  - c. Positive Control (fortified tissue)
  - d. Samples
  - e. Positive Control (fortified tissue)
7. Minimum Level of Applicability (MLA): 15 ppb

**J. WORKSHEET**

{Reserved}

**K. APPENDIX**


1. Chromatograms
  - a. Negative control



**United States Department of Agriculture  
Food Safety and Inspection Service, Office of Public Health Science**

CLG-CBX4.00	Page 12 of 12	
Title: Screen of Carbadox Metabolite QCA by HPLC-MS-MS		
Revision:NA	Replaces: NA	Effective: 03/28/2011

b. Positive control at 15 ppb



2. Reference:

Hutchinson, M.J., Young, P.Y., Hewitt, S.A., Faulkner, D., Kennedy, D.G., "Development and validation of an improved method of the carbadox metabolite, quinoxaline- 2- carboxylic acid, in porcine liver using LC-electrospray MS-MS according to revised EU criteria for veterinary drug residue analysis", *The Analyst* 2002, 127, 342-346.

**L. APPROVALS AND AUTHORITIES**

1. Approvals on file.
2. Issuing Authority: Director, Laboratory Quality Assurance Division.