

**United States Department of Agriculture
Food Safety and Inspection Service, Office of Public Health Science**

CLG-PBZ3.02		Page 1 of 13
Title: Screening of Phenylbutazone in Bovine Kidney by ELISA		
Revision: 02	Replaces: CLG-PBZ3.01	Effective: 08/18/2008

Contents

A.	INTRODUCTION	2
B.	EQUIPMENT	2
C.	REAGENTS AND SOLUTIONS	3
D.	STANDARDS.....	4
E.	SAMPLE PREPARATION.....	5
F.	ANALYTICAL PROCEDURE	5
G.	CALCULATIONS	7
H.	SAFETY INFORMATION AND PRECAUTIONS	8
I.	QUALITY ASSURANCE PLAN	9
J.	WORKSHEET.....	11
K.	APPENDIX.....	13
L.	APPROVALS AND AUTHORITIES.....	13

**United States Department of Agriculture
Food Safety and Inspection Service, Office of Public Health Science**

CLG-PBZ3.02	Page 2 of 13	
Title: Screening of Phenylbutazone in Bovine Kidney by ELISA		
Revision: 02	Replaces: CLG-PBZ3.01	Effective: 08/18/2008

A. INTRODUCTION

1. Theory

This procedure is used to screen bovine kidney samples for phenylbutazone residues at concentrations of 50 ppb or higher. Blended kidney tissues are extracted with a dilute acidic phosphate buffer and analyzed using a commercial Phenylbutazone ELISA kit (PBZ ELISA kit, Neogen Corporation). Samples exhibiting absorbances lower than a calculated decision level are identified as positives.

2. Applicability

The method is applicable to the detection of phenylbutazone in bovine kidney at levels ≥ 50 ppb.

B. EQUIPMENT

1. Apparatus

Equivalent apparatus may be substituted for the following:

- a. Balance, analytical - Model No. MT5, sensitive to 0.0001 g, Mettler.
- b. Balance, top loader - Model No. PM2000, sensitive to 0.01g, Mettler.
- c. Vortex mixer - Super Mixer Model No.1290, Lab-Line.
- d. Shaker, platform (two speed) - Eberhard.
- e. Volumetric dispenser, bottle top - Capable of delivering from 1 to 10 mL, Brinkmann.
- f. Stirrer/hot plate - Model No PC-4420, Corning.
- g. Centrifuge - International Equipment Company B-22M Superspeed Refrigerated Centrifuge with Rotor No. 876, International Equipment Company, Clinical Centrifuge, Cat. No. 20671-007, VWR Scientific.
- h. Eppendorf pipettors - Variable volume pipettes, 2-20 μ L (Cat. No. 05-402-46), 10-100 μ L (Cat. No. 05-402-48), 50-200 μ L (Cat. No. 05-402-49), 100-1000 μ L (Cat. No. 05-402-50) and 500-5000 μ L (Cat. No. 022472151),.
- i. epMotion 5075 - Liquid Handling Workstation from Eppendorf (Cat. No. 960020006), TS 50 – 1 to 50 μ L dispensing tool (Cat. No. 960001010), TS 1000 – 40 to 1000 μ L dispensing tool (Cat. No. 960001036), TM 1000-8 – 40 to 1000 μ L 8-channel dispensing tool (Cat. No. 960001061).
- j. pH meter - Orion 601A Calibrated at pH 4 and 7, readable to 0.01 pH.
- k. Transferpette, Multichannel Pipettes - Brinkmann, 50-300 μ L #50-08-030-7 (Cat.

**United States Department of Agriculture
Food Safety and Inspection Service, Office of Public Health Science**

CLG-PBZ3.02	Page 3 of 13	
Title: Screening of Phenylbutazone in Bovine Kidney by ELISA		
Revision: 02	Replaces: CLG-PBZ3.01	Effective: 08/18/2008

No. 53512-376, VWR).

- i. Glassware - Volumetric glassware includes 50-mL, 100-mL, 1000-mL, and 2000-mL flasks. CLASS A.
 - m. Pipettes - 25 mL glass pipettes. CLASS A.
 - n. Test tubes (for use with Liquid Handling Workstation) - 12mm x 75mm disposable Borosilicate glass culture tubes, (Kimble).
 - o. Centrifuge tubes - 50-mL polyallomer tube with polypropylene screw closure, Cat. No. 3139-0050, Nalge Company.
2. Instrumentation
- a. Plate Reader: Biotek Autoreader ELx 808 - Equipped with 650 nm filter and a printer (ELISA Technologies).

C. REAGENTS AND SOLUTIONS

1. Reagents

Note: Equivalent reagents may be substituted for the following:

- a. Test kit – Phenylbutazone Enhanced ELISA Kit, Cat. No. 104710-I or 104715-I (bulk kit) (Neogen Corporation, ELISA Technologies Division).

Contents of the bulk kit include the following:

- i. EIA Buffer (200 mL).
 - ii. 10x Wash Buffer Concentrate (100 mL). Dilute before use according to manufacturer's instructions.
 - iii. K-Blue Substrate (100 mL). Stabilized 3,3',5,5'-Tetramethylbenzidine (TMB) plus hydrogen peroxide in a single bottle. Light sensitive.
 - iv. Drug-Enzyme Conjugate (1mL). Drug-horseradish peroxidase concentrate. Dilute before use according to manufacturer's instructions.
 - v. Precoated Plates (5). Each Costar plate consists of 96 wells in strips of 8 breakaway wells coated with anti-drug antisera. The plates are ready for use. Do not wash until the sample/drug-conjugate incubation is complete.
- Note: The test kits must be stored in a refrigerator at 2 - 8 °C. Do not use past expiration date.
- b. Sodium Phosphate, dibasic, anhydrous. Na₂HPO₄, Cat. No. 3828-01, J.T.Baker.
 - c. Potassium Phosphate, Monobasic crystals. KH₂PO₄, Cat. No. 3246-01, J.T.Baker.
 - d. Deionized water.

**United States Department of Agriculture
Food Safety and Inspection Service, Office of Public Health Science**

CLG-PBZ3.02	Page 4 of 13	
Title: Screening of Phenylbutazone in Bovine Kidney by ELISA		
Revision: 02	Replaces: CLG-PBZ3.01	Effective: 08/18/2008

- e. Methanol, LC-Grade, 4L, Cat. No. H488-10, Mallinckrodt.
- f. 1N HCl, 1L, Cat No. SA48-1, Fisher Scientific.
- g. 1N NaOH, 1L, Cat. No. SS266-1, Fisher Scientific.
- h. Buffer solution, pH 4.01, Cat. No. 238217, Hamilton.
- i. Buffer solution, pH 7.00, Cat. No. 238218, Hamilton.

2. Solutions

Note: Equivalent solutions may be substituted for the following:

- a. Extraction Buffer (0.2M Phosphate)

Weigh 13.6 g of potassium phosphate monobasic, and 14.2 g of sodium phosphate dibasic into a 1 L class A volumetric flask. Add approximately 950 mL of deionized water to dissolve reagents. Adjust pH to 6.8 ± 0.1 with 1N HCl or 1N NaOH solutions. Dilute to final volume with deionized water. Prepare fresh extraction buffer when it becomes cloudy. The pH should be checked periodically to verify that it is 6.8 ± 0.1 .
- b. Horseradish peroxidase (HRP) enzyme conjugate (1:180)– EIA buffer dilution (1 + 179). Prepare in a ratio of 1 μ L HRP conjugate to 179 μ L EIA buffer each day of use or per manufacturer's instructions.
- c. Wash buffer, dilution – Dilute 1:9 with deionized water or per manufacturer's instructions. This solution is assigned the same expiration date as the concentrated wash buffer and shall be stored at 2 - 8 °C.

D. STANDARDS

1. Source

Phenylbutazone standard is available from MP Biomedicals. (Cat. No. 153567). Other sources for the standard material may be acceptable.

2. Preparation

Note: Equivalent solutions may be substituted for the following:

- a. Stock PBZ Standard Solution, 500 μ g/mL in methanol:

Weigh 50.0 ± 0.1 mg PBZ standard into a 100-mL volumetric flask. Dissolve and bring to volume with methanol.
- b. Intermediate PBZ Standard Solution, 250 μ g/mL in methanol:

**United States Department of Agriculture
Food Safety and Inspection Service, Office of Public Health Science**

CLG-PBZ3.02	Page 5 of 13	
Title: Screening of Phenylbutazone in Bovine Kidney by ELISA		
Revision: 02	Replaces: CLG-PBZ3.01	Effective: 08/18/2008

Pipet 25 mL PBZ stock standard (D.2.a.) into a 50-mL volumetric flask and bring to volume with methanol.

c. Working PBZ Standard Solution, 300 ng/mL:

Dilute 60 μ L of intermediate standard (D.2.b.) to 50 mL with extraction buffer (C.2.a.) in a 50-mL volumetric flask.

Note: The working standard solution is used to generate the standard curve and for sample fortification.

d. External Standard Solutions (0, 15, and 30 ng/mL):

i. 0 ng/mL (0 ppb): extraction buffer solution (C.2.a).

ii. 15 ng/mL (50 ppb): Dilute 50 μ L of Working Standard (D.2.c.) to 1 mL with extraction buffer solution (C.2.a).

iii. 30 ng/mL (100 ppb): Dilute 100 μ L of Working Standard (D.2.c.) to 1 mL with extraction buffer solution (C.2.a).

3. Storage and Stability

Store Standards D.2.a. and D.2.b. at 2 - 8 °C and Standard D.2.c. at room temperature. Stock and intermediate standards are stable for 2 months. The working PBZ standard solution (D.2.c) and External Standard Solutions (D.2.d.ii and D.2.d.iii) must be prepared fresh daily.

E. SAMPLE PREPARATION

1. Kidney Samples

Intact kidney samples (one/cow), received cool or frozen, are halved, with one half of the kidney blended and the other half stored frozen, whenever sufficient sample is available.

2. Negative control tissue

Negative control tissue is prepared by forming a composite of at least six different bovine kidney samples which were previously found to be screened negative. The tissues are combined and blended to ensure homogeneity. Refer to I.5.c.i. for storage conditions.

F. ANALYTICAL PROCEDURE

1. Extraction Procedure

a. Preparation of blank samples for establishing the decision level

i. Weigh 6 known kidney blanks (4 ± 0.1 g each) into 50 mL polyallomar centrifuge tubes.

**United States Department of Agriculture
Food Safety and Inspection Service, Office of Public Health Science**

CLG-PBZ3.02	Page 6 of 13	
Title: Screening of Phenylbutazone in Bovine Kidney by ELISA		
Revision: 02	Replaces: CLG-PBZ3.01	Effective: 08/18/2008

- ii. Add 2.0 mL of extraction buffer (C.2.a.) to each of six tubes.
- iii. Continue the process following steps F.1.c.iii. through F.1.c.v.
- b. Preparation of 50 ppb control sample and negative control
 - i. Weigh one known kidney blank (4 ± 0.1 g each) into two 50-mL polyallomar centrifuge tubes.
 - ii. For the 50 ppb control sample, add 1.0 mL of 300 ng/mL working standard (D.2.c.) plus 1.0 mL of extraction buffer (C.2.a.) to the tube. For the 0 ppb negative control, add 2.0 mL of extraction buffer (C.2.a.) to the tube.

Note: The fortification level is based on the concentration of PBZ added to the tissue divided by the total g of product (4 g of tissue plus 2 g of extraction buffer or standard solution = 6 g total). Application of 20 μ L of these solutions on each well represents 0 ng and 3 ng of PBZ, respectively, for 0 and 50 ppb fortifications.
 - iii. Continue the process following steps F.1.c.iii. through F.1.c.v.
- c. Extraction of samples
 - i. Weigh 4 ± 0.1 g of blended sample into a 50 mL polyallomar centrifuge tube.
 - ii. Add 2.0 mL of extraction buffer (C.2.a.) to each tube.
 - iii. Vortex vigorously for at least 1 minute and sonicate for 5 minutes.
 - iv. Centrifuge tubes for 15 min at ~15,000 RPM at 4 °C.
 - v. Extracts are ready for plating. Keep extracts refrigerated until beginning step F.2.

2. ELISA

Note: The following steps are based on kit manufacturer instructions and may be subject to change. If any discrepancies exist, follow the current manufacturer instructions.

- a. Apply in duplicate 20 μ L extracts of each external standard(s), control(s), and sample(s) into individual ELISA plate wells.

Note: The six blank extracts must be distributed randomly in duplicate into 12 wells covering the area of the plate used in order to facilitate estimation of a decision level (see Section G).

- b. Add 180 μ L HRP: EIA solution (C.2.b) to each well. Mix the solutions by gently vibrating the plate on a flat surface.
- c. Cover plate to avoid possible dust/dirt contamination.

**United States Department of Agriculture
Food Safety and Inspection Service, Office of Public Health Science**

CLG-PBZ3.02	Page 7 of 13	
Title: Screening of Phenylbutazone in Bovine Kidney by ELISA		
Revision: 02	Replaces: CLG-PBZ3.01	Effective: 08/18/2008

- d. Incubate for one hour at room temperature. Shake gently at least twice during the incubation period.
- e. Invert the plate after the incubation period, to remove matrix solutions.
- f. Wash the wells 4 times using 300 μ L/well per wash of diluted wash buffer (C.2.c). Tap the inverted plate on a paper towel between washings and make sure that the plate is free from liquid and bubbles prior to adding K-Blue substrate.
- g. Add 150 μ L K-Blue substrate (another reagent supplied in the test kit) to each well. Allow the reaction to proceed for 15–30 minutes with intermittent gentle shaking of plate, especially before taking an optical density (absorbance) reading.

3. Instrumental Settings

Evaluate the plate using a plate reader at 650 nm.

Note: Optimum total incubation time will be that where the absorbance reading for the plate blank just begins to plateau. This may be established by recording plate absorbances at 5 minute intervals until the absorbance reading of the 0 external standard no longer shows a significant increase.

G. CALCULATIONS

Evaluate sample results based on absorbance values for unfortified control tissues. The extent of the color development in each well is inversely proportional to the amount of drug in the sample or control.

Calculate the mean and standard deviation (SD) for the absorbance readings of the six unfortified control tissue replicates (see F.2.a.). Use these to calculate a decision level (DL) using the formula $DL = \text{Mean} - 2.78 \times SD$.

A sample run using duplicate aliquots is screen positive if it meets one of the following conditions:

- a. Both absorbance readings are less than the decision level.
- b. One reading is less than the decision level and the average of both readings is less than the decision level and the relative percent difference between readings is less than 10%. Determined as follows:

Calculate the absolute difference between the two measurements and divide by two. Multiply this number by 100. Divide this product by the mean of the two measurements (resultant calculation is expressed as percentage).

Note: If one well has an absorbance greater than the decision level and the relative difference between duplicate wells is greater than 10%, the sample must be analyzed again.

**United States Department of Agriculture
Food Safety and Inspection Service, Office of Public Health Science**

CLG-PBZ3.02	Page 8 of 13	
Title: Screening of Phenylbutazone in Bovine Kidney by ELISA		
Revision: 02	Replaces: CLG-PBZ3.01	Effective: 08/18/2008

H. SAFETY INFORMATION AND PRECAUTIONS

1. Required Protective Equipment — Lab coat, safety glasses and gloves.
2. Hazards

<i>Procedure Step</i>	<i>Hazard</i>	<i>Recommended Safe Procedures</i>
Methanol	Flammable; may produce toxic effects to skin, eyes and respiratory system.	Use reagents in an efficient fume hood away from all electrical devices and open flames.
Acid and base	Corrosive	Wear gloves and safety glasses

3. Disposal Procedures

<i>Procedure Step</i>	<i>Hazard</i>	<i>Recommended Safe Procedures</i>
Methanol	See above	Collect waste in tightly sealed container and store away from non-compatibles in a cool, well ventilated, flammable liquid storage area/cabinet for disposal in accordance with local, state, and Federal regulations.
Acid and base	Can cause burns on skin and eye injury	Neutralize the acid or base for disposal down the drain in accordance with local, state, and Federal regulations.
Used plates	None	Dispose of plates in accordance with local, state, and Federal regulations.

**United States Department of Agriculture
Food Safety and Inspection Service, Office of Public Health Science**

CLG-PBZ3.02	Page 9 of 13	
Title: Screening of Phenylbutazone in Bovine Kidney by ELISA		
Revision: 02	Replaces: CLG-PBZ3.01	Effective: 08/18/2008

I. QUALITY ASSURANCE PLAN

1. Performance Standard

<i>Analyte</i>	<i>Analytical Range</i>	<i>Acceptable Recovery</i>	<i>Acceptable Repeatability (CV)</i>
PBZ	≥ 50 ppb	NA	NA

- a. Each sample set must meet all the following criteria before data can be reported:
 - i. The negative control tissue absorbances should be between 1.0 and 2.0 absorbance units vs. air.
 - ii. The external standard curve absorbance values continuously increase from the 100 ppb to the 0 ppb concentration.
 - iii. The 0 ppb standard and 0 ppb fortified control wells must screen negative.
 - iv. The CV calculated for the six negative control tissue replicates (12 measurements) used to establish the decision level must be ≤ 20%.
- b. No false negatives at ≥ 50 ppb fortified tissue.

2. Critical Control Points and Specifications

Record

Acceptable Control

- | | |
|---|---|
| <ol style="list-style-type: none"> a. After washing the ELISA plate make sure that the plate is free of bubbles and liquid prior to proceeding to next step of adding K-blue solution (F.2.g). | <ol style="list-style-type: none"> Evacuate plate of all liquid or bubbles prior to the addition of the K-blue solution. |
|---|---|

3. Readiness To Perform

- a. Familiarization
 - i. Phase I: Standards-Duplicate standard curve (0, 15, and 30 ng/mL) on each of 3 different days.
 - ii. Phase II: Fortified samples-3 sets of 10 blank and 10 fortified samples at 50 ppb level over a period of 3 different days.
Note: Phase I and Phase II may be performed concurrently.
 - iii. Phase III: Check samples for analyst accreditation.
 - (i) The analyst is to receive 30 bovine kidney samples. The sample fortifications, including the number of blanks, are to be blind to the

**United States Department of Agriculture
Food Safety and Inspection Service, Office of Public Health Science**

CLG-PBZ3.02	Page 10 of 13	
Title: Screening of Phenylbutazone in Bovine Kidney by ELISA		
Revision: 02	Replaces: CLG-PBZ3.01	Effective: 08/18/2008

analyst. At least 10 of the 30 samples should be blank and the rest fortified at the 50 ppb level. These samples will be prepared and provided by Supervisor or QA Manager (QAM). All samples should be analyzed in duplicate wells. The samples must be randomized throughout the set.

- (ii) An external curve at 0 ng/mL, 15 ng/mL and 30 ng/mL (equivalent to 0, 50, and 100 ppb, respectively) must be run to help monitor plate acceptability.
- (iii) Report analytical findings to the Laboratory Quality Assurance Manager.
- (iv) Authorization from the QAM and Supervisor is required to commence official sample analysis.

b. Acceptability criteria.

Refer to I. 1. a. and b.

4. Intralaboratory Check Samples

a. System, minimum contents.

- i. Frequency: At least 1 weekly per analyst when samples are analyzed.
- ii. Records of check sample results are to be maintained.

b. Acceptability criteria.

Refer to I. 1.

If unacceptable values are obtained, then:

- i. Stop all official analyses by that analyst.
- ii. Take corrective action.

5. Sample Acceptability and Stability

- a. Matrix: Bovine kidney.
- b. Sample receipt size: One intact kidney.
- c. Sample storage:
 - i. Condition: Frozen (< -10 °C).

6. Sample Set

- a. Each sample set contains:
 - i. Six different control tissue blanks to establish the decision level.
 - ii. One negative control.

**United States Department of Agriculture
Food Safety and Inspection Service, Office of Public Health Science**

CLG-PBZ3.02	Page 11 of 13	
Title: Screening of Phenylbutazone in Bovine Kidney by ELISA		
Revision: 02	Replaces: CLG-PBZ3.01	Effective: 08/18/2008

- iii. Fortified tissue at 50 ppb.
- iv. Standards at 0, 15, and 30 ng/mL.
- v. Samples.

Note: At least 5 % of the sample set must be comprised of QC samples (50 ppb fortified(s) and/or Check Sample(s)).

- 7. Sensitivity
 - a. Minimum proficiency level (MPL): 50 ppb.

J. WORKSHEET

Following is an example worksheet.

**United States Department of Agriculture
Food Safety and Inspection Service, Office of Public Health Science**

CLG-PBZ3.02	Page 12 of 13	
Title: Screening of Phenylbutazone in Bovine Kidney by ELISA		
Revision: 02	Replaces: CLG-PBZ3.01	Effective: 08/18/2008

CLG-PBZ3 PHENYL BUTAZONE SCREENING FORM

Weighed (Init./AC/Date): _____

Test Portion Holding FRZ: _____ Init./AC/Date: _____

Extracted (Init./AC/Date): _____

Reader Analysis (Init./AC/Date): _____

Set Number: _____

Method: CLG-PBZ3.XX

Reviewed By Initials/AC/Date: _____

Standards/Reagents Used	MWL ID #
Fortification standard (s)	
External standard (s)	
0.2M phosphate extraction buffer	
PBZ ELISA kit	
1:180 HRP-EIA diluted solution	
Diluted wash buffer	
K-Blue substrate	
Equipment #	MWL ID #
Refrigerator (REF)	
Freezer (FRZ)	
Balance (BAL)	
Micropipettors (MCP)	
Sonicator (WAT)	
Shaker (TIM)	
Centrifuge (CNT)	
Optical density reader (SPC)	

Position Number	Lab Number	Form Number	Sample Wt. 4.00 ± 0.10 g	Results (+) or (-)
1				
2				
3				
4				
5				
6				
7				
8				
9				
10				
11				
12				
13				
14				
15				
16				
17				
18				
19				
20				
21				
22				
23				
24				

Position Number	Lab Number	Form Number	Sample Wt. 4.00 ± 0.10 g	Results (+) or (-)
25				
26				
27				
28				
29				
30				
31				
32				
33				
34				
35				
36				
37				
38				
39				
40				
41				
42				
43				
44				
45				
46				
47				
48				

**United States Department of Agriculture
Food Safety and Inspection Service, Office of Public Health Science**

CLG-PBZ3.02	Page 13 of 13	
Title: Screening of Phenylbutazone in Bovine Kidney by ELISA		
Revision: 02	Replaces: CLG-PBZ3.01	Effective: 08/18/2008

K. APPENDIX

1. Susan B. Clark, Sherri B. Turnipseed, Gene J. Nandrea, Mark R. Madson, Emma R. Singleton, Jeffrey A. Hurlbut, John N. Sofos, and Craig E. Shultz. "Identification and Confirmation of Flunixin Meglumine and Phenylbutazone Residues in Animal Kidney by ELISA Screening and Liquid Chromatography Mass Spectrometry." LIB 4246, Food and Drug Administration, Denver, CO. Vol. 17, No. 5, May 2001.

L. APPROVALS AND AUTHORITIES

Approvals on file.

Issuing Authority: Director, Laboratory Quality Assurance Division.