

**United States Department of Agriculture  
Food Safety and Inspection Service, Office of Public Health Science**

|                                                                                       |                       |                     |
|---------------------------------------------------------------------------------------|-----------------------|---------------------|
| CLG-FLQ2.01                                                                           |                       | Page 1 of 20        |
| Title: Confirmation of Fluoroquinolone Antibiotics by HPLC Ion Trap Mass Spectrometry |                       |                     |
| Revision: 01                                                                          | Replaces: CLG-FLQ2.00 | Effective: 04/26/21 |

**Contents**

|    |                                          |    |
|----|------------------------------------------|----|
| A. | INTRODUCTION .....                       | 2  |
| B. | EQUIPMENT .....                          | 2  |
| C. | REAGENTS AND SOLUTIONS .....             | 3  |
| D. | STANDARD(S) .....                        | 5  |
| E. | SAMPLE RECEIPT AND PREPARATION .....     | 6  |
| F. | ANALYTICAL PROCEDURE .....               | 6  |
| G. | DECISION CRITERIA / CALCULATIONS .....   | 10 |
| H. | SAFETY INFORMATION AND PRECAUTIONS ..... | 11 |
| I. | QUALITY ASSURANCE PLAN .....             | 11 |
| J. | APPENDIX.....                            | 12 |
| K. | APPROVALS AND AUTHORITIES.....           | 20 |

**United States Department of Agriculture  
Food Safety and Inspection Service, Office of Public Health Science**

|                                                                                       |                       |                     |
|---------------------------------------------------------------------------------------|-----------------------|---------------------|
| CLG-FLQ2.01                                                                           | Page 2 of 20          |                     |
| Title: Confirmation of Fluoroquinolone Antibiotics by HPLC Ion Trap Mass Spectrometry |                       |                     |
| Revision: 01                                                                          | Replaces: CLG-FLQ2.00 | Effective: 04/26/21 |

**A. INTRODUCTION**

1. Background / Summary of Procedure

Fluoroquinolone antibiotics (FLQs) are extracted from homogenized tissue using a liquid/liquid technique. Upon concentration of the extracts, eight fluoroquinolones are analyzed by HPLC ion trap mass spectrometry. Two fluoroquinolones (desethylene ciprofloxacin and desmethyl danofloxacin) are analyzed using HPLC/MS<sup>2</sup>, whereas six other fluoroquinolones (difloxacin, enrofloxacin, norfloxacin, danofloxacin, ciprofloxacin, sarafloxacin) are analyzed by HPLC/MS<sup>3</sup>.

2. Applicability

This method is suitable for the confirmation of desethylene ciprofloxacin, difloxacin, enrofloxacin, norfloxacin, danofloxacin, desmethyl danofloxacin, ciprofloxacin and sarafloxacin in bovine liver and muscle at levels  $\geq 25$  ppb.

*Note: Refer to 21CFR for tolerance values set by FDA and 40CFR for tolerance values set by EPA.*

*Note: This method may be performed using standards/solutions that contain fewer analytes than the method is applicable for, if the excluded analytes will not be included in the reported results.*

**B. EQUIPMENT**

*Note: Equivalent equipment may be substituted.*

1. Apparatus

- a. Waring Blender - Model BLH 120, Waring Inc.
- b. Robot Coupe® - Model RSI 3Y-1, Robot Coupe Inc.
- c. Centrifuge - Model KR22i, Jouan.
- d. Homogenizer - Model IKA, Ultra Turrax.
- e. Vortex mixer - Fisher-Genie 2.
- f. Pasteur pipettes - borosilicate glass, 5.75 inches.
- g. Nylon syringe filters - 13 mm, 0.22 $\mu$ m disposable, Cat. No. 9445622, Xpertext.
- h. Balance - accurate to 0.0001 g, Cat. No. MT5, Mettler.
- i. Nitrogen evaporator - Turbovap LV, Zymark.
- j. Plasticware - 50 mL polypropylene centrifuge tubes, 30 x115 mm, Cat. No. 352070, Falcon.
- k. Autosampler vials - 750  $\mu$ L, plastic, Cat. No. 951501, Xpertext.

**United States Department of Agriculture  
Food Safety and Inspection Service, Office of Public Health Science**

|                                                                                       |                       |                     |
|---------------------------------------------------------------------------------------|-----------------------|---------------------|
| CLG-FLQ2.01                                                                           | Page 3 of 20          |                     |
| Title: Confirmation of Fluoroquinolone Antibiotics by HPLC Ion Trap Mass Spectrometry |                       |                     |
| Revision: 01                                                                          | Replaces: CLG-FLQ2.00 | Effective: 04/26/21 |

- l. Micropipetters - covering the range from 10  $\mu$ L – 5000  $\mu$ L.
  - m. Vibrax mixer - Model VX8, Janke & Kunkel.
  - n. Balance - accurate to 0.01 g, Cat. No. PB302, Mettler.
2. Instrumentation
- a. Ion trap mass spectrometer - Thermo LTQ XL Linear Ion Trap with ESI interface and Windows NT ver.4.0- LCQ Xcalibur data system, or equivalent.
  - b. LC system - Quaternary pump equipped with degassing capability and autosampler.
  - c. LC column - Phenyl 3 x 150 mm containing 3.5  $\mu$ m particles, Cat. No. XDB, Zorbax.
  - d. Guard column - 2.1 mm x 12.5 mm containing 5  $\mu$ m particles, Eclipse XDB-C8.

**C. REAGENTS AND SOLUTIONS**

*Note: Equivalent reagents / solutions may be substituted.*

1. Reagents
- a. Ethyl ether - ACS Grade, 99%, Aldrich.
  - b. Water, LC grade - House distilled water passed through Waters MilliQ deionization system.
  - c. Acetonitrile - UV Grade, Cat. No. 015-4, Burdick & Jackson.
  - d. Sodium chloride - ACS Grade, Cat. No. 3624-01, Baker.
  - e. Ammonium hydroxide 28% - Cat. No. 38,053-9, Aldrich.
  - f. Hexane - Omnisolv. Cat. No. HX0296-1, EM.
  - g. Sodium phosphate dibasic, heptahydrate - ACS Grade, Cat. No. S-9390, Sigma.
  - h. Formic acid - ACS Grade, Cat. No. F-4636, Sigma.
  - i. Sodium phosphate monobasic, monohydrate - ACS Grade, Cat. No. S-9638, Sigma.
  - j. Sodium hydroxide - ACS Grade.
  - k. Methanol - HPLC grade, Mallinckrodt.
2. Solutions
- a. 1M sodium chloride solution:  
Transfer 58.45 g NaCl to a 1000 mL volumetric flask. Dissolve and dilute to volume with water.

**United States Department of Agriculture  
Food Safety and Inspection Service, Office of Public Health Science**

|                                                                                       |                       |                     |
|---------------------------------------------------------------------------------------|-----------------------|---------------------|
| CLG-FLQ2.01                                                                           | Page 4 of 20          |                     |
| Title: Confirmation of Fluoroquinolone Antibiotics by HPLC Ion Trap Mass Spectrometry |                       |                     |
| Revision: 01                                                                          | Replaces: CLG-FLQ2.00 | Effective: 04/26/21 |

- b. 0.2M monobasic sodium phosphate monohydrate:  
Transfer 27.6 g  $\text{NaH}_2\text{PO}_4$  to a 1000 mL volumetric flask. Dissolve and dilute to volume with water. Solution is stable for 6 months at 2-8 °C unless the note above applies.
- c. 0.2M dibasic sodium phosphate heptahydrate:  
Transfer 53.65 g  $\text{Na}_2\text{HPO}_4$  to a 1000 mL volumetric flask. Dissolve and dilute to volume with water. Solution is stable for 6 months at 2-8 °C unless the note above applies.
- d. 30% sodium hydroxide:  
Transfer 60 g NaOH to a 200 mL volumetric flask and dilute to volume with water. Solution is stable for 6 months at 2-8 °C unless the note above applies.
- e. 0.1M phosphate buffer, pH 9.0 (Buffer A):  
Transfer 27 mL of 0.2M monobasic sodium phosphate monohydrate (a) and 473 mL of 0.2M dibasic sodium phosphate heptahydrate (b) to a 1000 mL beaker. Add water to approx. 900 mL. Adjust to pH 9.0 with 30% NaOH (c). Transfer to a 1000 mL volumetric flask and dilute to volume with water. Solution is stable for 6 months at 2-8 °C unless the note above applies.
- f. 0.03M sodium hydroxide:  
Transfer 1.2 g NaOH to a 1000 mL volumetric flask and dilute to volume with water.
- g. Mobile phase A (15/85 acetonitrile/water containing 1% formic acid):  
To a 1000 mL graduated cylinder add 150 mL acetonitrile and 840 mL water. Mix and add 10 mL formic acid. Prepare fresh water from water purification system at the time of each run. Mix well.
- h. Mobile phase B (20/80 acetonitrile/water containing 1% formic acid):  
To a 1000 mL graduated cylinder add 200 mL acetonitrile and 790 mL water. Mix and add 10 mL formic acid. Mix well.
- i. Mobile phase C (80/20 acetonitrile/water containing 1% formic acid):  
To a 1000 mL graduated cylinder add 800 mL acetonitrile and 190 mL water. Mix and add 10 mL formic acid. Mix well.
- j. Mobile phase D (60/40 acetonitrile/water):  
To a 1000 mL graduated cylinder add 600 mL acetonitrile and 400 mL water. Mix well.

**United States Department of Agriculture  
Food Safety and Inspection Service, Office of Public Health Science**

|                                                                                       |                       |                     |
|---------------------------------------------------------------------------------------|-----------------------|---------------------|
| CLG-FLQ2.01                                                                           | Page 5 of 20          |                     |
| Title: Confirmation of Fluoroquinolone Antibiotics by HPLC Ion Trap Mass Spectrometry |                       |                     |
| Revision: 01                                                                          | Replaces: CLG-FLQ2.00 | Effective: 04/26/21 |

**D. STANDARD(S)**

*Note: Equivalent standards / solutions may be substituted. Purity and counterions are to be taken into account when calculating standard concentrations. In-house prepared standards shall be assigned an expiration date that is no later than the stability stated in the method.*

1. Standard Information

- a. Desethylene ciprofloxacin: Bayer Corp., West Haven, CT.
- b. Difloxacin HCl: Abbott Labs, North Chicago, IL.
- c. Enrofloxacin: Bayer Corp., West Haven, CT.
- d. Norfloxacin: Sigma-Aldrich, St. Louis, MO.
- e. Ciprofloxacin: United States Pharmacopoeia (USP).
- f. Desmethyl danofloxacin: Pfizer Pharmaceuticals, Groton, CT.
- g. Sarafloxacin HCl: Abbott Labs, North Chicago, IL.
- h. Danofloxacin Mesylate: Pfizer Pharmaceuticals, Groton, CT.

2. Preparation of Standard Solution(s)

- a. Individual drug stock standard solutions (100 µg/mL):  
Using vendor's stated purity, or water and salt content, calculate the amount of material which contains 5 mg drug. Weigh out approximately this amount, accurately recording weight to nearest 0.1 mg. Transfer to 50 mL glass volumetric flask and dilute to mark with 0.03M NaOH (C.2.f). Calculate exact concentration based on purity and actual weight. Standard is stable for 6 months at 2-8 °C unless the note above applies.
- b. Mixed working standard solution (2 µg/mL):  
Add 1.0 mL of each of the above drug stock solutions to a 50 mL volumetric flask and dilute to mark with Buffer A (C.2.e.). Standard is stable for 1 month at 2-8 °C unless the note above applies.
- c. External standard solution (12.5 ng/mL):  
Transfer 12.5 µL of the mixed working standard solution (b) to a 50 mL plastic centrifuge tube and add 2.0 mL of Buffer A (C.2.e). Prepare daily.
- d. System suitability standard (12.5 ng/mL):  
Transfer 12.5 µL of the mixed working standard solution (b) to a 50 mL plastic centrifuge tube and add 2.0 mL of methanol. Prepare daily.

**United States Department of Agriculture  
Food Safety and Inspection Service, Office of Public Health Science**

|                                                                                       |                       |                     |
|---------------------------------------------------------------------------------------|-----------------------|---------------------|
| CLG-FLQ2.01                                                                           | Page 6 of 20          |                     |
| Title: Confirmation of Fluoroquinolone Antibiotics by HPLC Ion Trap Mass Spectrometry |                       |                     |
| Revision: 01                                                                          | Replaces: CLG-FLQ2.00 | Effective: 04/26/21 |

**E. SAMPLE RECEIPT AND PREPARATION**

Freshly collected samples must be kept cold before and during shipping to laboratory. Once received at laboratory, samples must be frozen (< -10 °C) prior to processing if they cannot be prepared on the day of receipt. If sample is frozen, allow to thaw, but keep as cold as possible. Dissect away fat and connective tissue from liver and muscle. Homogenize liver in a Waring blender and muscle in a Robot Coupe®.

Note: After each homogenization, rinse the blending jars with tap water and dry.

**F. ANALYTICAL PROCEDURE**

**1. Preparation of Controls and Samples**

- a. Weigh  $1.0 \pm 0.10$  g of homogenized tissue samples into a 50 mL disposable polypropylene centrifuge tube.
- b. Weigh two 1.0 g portions of blank tissue into 50 mL polypropylene centrifuge tubes, allow tissue to thaw and do the following described below:
  - i. Prepare one each for a negative control and a positive control.
  - ii. Prepare a 25 ppb positive control, fortify with 12.5  $\mu$ L of the 2  $\mu$ g/mL mixed drug solution (D.2.b) and vortex 30 min on a Vibrax mixer. Store in a dark cabinet (room temp.) for 0.5 hr to allow time for the FLQ's to interact with matrix. Weigh one additional portion for an intra-laboratory check sample if necessary.

**2. Extraction Procedure**

- a. To the samples and controls add approx. 3 mL acetonitrile and 0.25 mL conc. ammonium hydroxide to each tube.
- b. Homogenize tube contents for approx. 20 sec. Rinse homogenizer tip with approx. 0.5 mL of acetonitrile directly into the sample tube. Centrifuge at approximately 4000 rpm (2800 g.) for 10 min.
- c. Decant supernatant into another 50 mL polypropylene centrifuge tube.
- d. Add 0.75 mL of water to the pellet from the centrifugation and repeat steps a-b, ensure the pellet breaks by vortexing for 30 sec, centrifuge at approximately 4000 rpm for 10 min. Add the supernatant to the tube in step c.
- e. Add approx. 3 mL ethyl ether, 3 mL hexane and 0.25 mL 1M NaCl to the combined supernatants for each sample.
- f. Vortex the tubes for 15 sec.
- g. Using a disposable Pasteur pipette or aspirate, remove as much of the top layer as possible and discard. There may be 3 layers at this point but discard only the uppermost layer.

**United States Department of Agriculture  
Food Safety and Inspection Service, Office of Public Health Science**

|                                                                                       |                       |                     |
|---------------------------------------------------------------------------------------|-----------------------|---------------------|
| CLG-FLQ2.01                                                                           |                       | Page 7 of 20        |
| Title: Confirmation of Fluoroquinolone Antibiotics by HPLC Ion Trap Mass Spectrometry |                       |                     |
| Revision: 01                                                                          | Replaces: CLG-FLQ2.00 | Effective: 04/26/21 |

- h. Evaporate the combined organic solution to near dryness (approximately 200 µL) in a Turbovap maintained at approximately 50 °C.
- i. Add 2.0 mL Buffer A (C.2.e) to each tube and vortex for approx 10-20 sec. Filter contents through a 0.2 µm nylon syringe filter into an autosampler vial.

3. Instrumental Settings

*Note: The instrument parameters may be optimized to ensure system suitability.*

a. Instrument Operating Parameters – LC System:

Note: Typical values listed below. Flows and elution gradient may be optimized, if necessary, for best separation and response.

- i. Install and degas mobile phases and install column and guard cartridge per manufacturers' instructions. Flush HPLC column with 20 column volumes (35 mL) of methanol, water, and 60/40 acetonitrile/water (C.2.j) prior to further use of the column. Set initial composition to flow 15/85 acetonitrile/water containing 1% formic acid (C.2.g) at 500 µL/min.
- ii. Set up the HPLC to run the following gradient:

| Time in min. | Flow in µL/min. | Mobile Phase A* | Mobile Phase B* | Mobile Phase C* |
|--------------|-----------------|-----------------|-----------------|-----------------|
| 0.00         | 200             | 100%            | 0%              | 0%              |
| 10.00        | 200             | 100%            | 0%              | 0%              |
| 18.00        | 200             | 0%              | 100%            | 0%              |
| 20.00        | 200             | 0%              | 100%            | 0%              |
| 22.00        | 200             | 0%              | 0%              | 100%            |
| 24.00        | 200             | 0%              | 0%              | 100%            |
| 27.00        | 200             | 100%            | 0%              | 0%              |
| 30.00        | 200             | 100%            | 0%              | 0%              |

A\* = C.2.g, B\* = C.2.h, C\* = C.2.i.

- iii. Set injection volume to 20 µL.
- iv. Use a needle wash step with methanol or water.
- b. Instrument Operating Parameters – Mass Spectrometer
  - i. Calibrate the Thermo LTQ XL Linear Ion Trap with ESI Interface according to the manufacturer's specifications.
  - ii. Set Capillary Temp to 200 °C.

**United States Department of Agriculture  
Food Safety and Inspection Service, Office of Public Health Science**

|                                                                                       |                       |                     |
|---------------------------------------------------------------------------------------|-----------------------|---------------------|
| CLG-FLQ2.01                                                                           |                       | Page 8 of 20        |
| Title: Confirmation of Fluoroquinolone Antibiotics by HPLC Ion Trap Mass Spectrometry |                       |                     |
| Revision: 01                                                                          | Replaces: CLG-FLQ2.00 | Effective: 04/26/21 |

- iii. Operate in Pos mode.
- iv. Flow inject the system suitability standard through a 5 µL loop and obtain the MS1 precursor ion centroids. The following settings should result in optimal ion intensities:

|                            |          |
|----------------------------|----------|
| Capillary temperature      | 200 °C   |
| APCI vaporizer temperature | 470 °C   |
| Sheath gas flow            | 50       |
| Aux gas flow               | 5        |
| Capillary voltage          | 2 V      |
| Tube lens offset           | -5 V     |
| Micro scans                | 2        |
| Ion time                   | 100 msec |
| Source current             | 5.00 µA  |
| *Spray Current             | >26      |

Note: \*Do not proceed with the sequence injection until the spray current is greater than 26

- c. Procedure for Instrumental Analysis of Samples, Controls, and Standards
  - i. Turn on pump and set up mass spectrometer. Equilibrate column in mobile phase at 200 µL/min for at least 30 min.
  - ii. Flow inject the system suitability standard through a 5 µL loop and obtain the MS1 precursor centroids. Using the MS1 mass assignments previously obtained, flow inject a sample under MS/MS conditions for each analyte and obtain the MS2 precursor ion centroids and accompanying collision energies.
  - iii. Inject the external standard through the HPLC system and determine the retention times of each FLQ. Isolate each peak within a suitable window for acquisition. After setting in the window segments, monitor the following ion transitions:

| Fluoroquinolone  | Scan Program     | Trans.          | Amp | Q     | Time | IsoW |
|------------------|------------------|-----------------|-----|-------|------|------|
| Desethylene Cip. | 306→(80-310)     | MS <sup>2</sup> | 35% | 0.250 | 30   | 2.0  |
| Norfloxacin      | 320→276→(75-330) | MS <sup>2</sup> | 30% | 0.350 | 30   | 2.0  |
|                  |                  | MS <sup>3</sup> | 35% | 0.250 | 30   | 2.0  |



**United States Department of Agriculture  
Food Safety and Inspection Service, Office of Public Health Science**

|                                                                                       |                       |                     |
|---------------------------------------------------------------------------------------|-----------------------|---------------------|
| CLG-FLQ2.01                                                                           |                       | Page 9 of 20        |
| Title: Confirmation of Fluoroquinolone Antibiotics by HPLC Ion Trap Mass Spectrometry |                       |                     |
| Revision: 01                                                                          | Replaces: CLG-FLQ2.00 | Effective: 04/26/21 |

|                 |                  |                 |     |       |    |     |
|-----------------|------------------|-----------------|-----|-------|----|-----|
| Ciprofloxacin   | 332→288→(75-340) | MS <sup>2</sup> | 35% | 0.350 | 30 | 2.0 |
|                 |                  | MS <sup>3</sup> | 35% | 0.250 | 30 | 2.0 |
| Danofloxacin    | 358→314→(85-370) | MS <sup>2</sup> | 35% | 0.350 | 30 | 2.0 |
|                 |                  | MS <sup>3</sup> | 35% | 0.250 | 30 | 2.0 |
| Desmethyl-Dano. | 344→(90-350)     | MS <sup>2</sup> | 40% | 0.250 | 30 | 2.0 |
| Enrofloxacin    | 360→316→(85-370) | MS <sup>2</sup> | 40% | 0.350 | 30 | 2.0 |
|                 |                  | MS <sup>3</sup> | 35% | 0.250 | 30 | 2.0 |
| Sarafloxacin    | 386→342→(90-400) | MS <sup>2</sup> | 35% | 0.350 | 30 | 2.0 |
|                 |                  | MS <sup>3</sup> | 38% | 0.250 | 30 | 2.0 |
| Difloxacin      | 400→356→(95-410) | MS <sup>2</sup> | 38% | 0.350 | 30 | 2.0 |
|                 |                  | MS <sup>3</sup> | 35% | 0.250 | 30 | 2.0 |

- iv. Inject the recovered standard and verify retention time, divert valve switching time, and spectral comparison to the external standard.
- v. Inject the sample extracts. In order to control carryover, precede each sample analysis with a blank buffer injection as needed.
- vi. As a test of retention time and instrument response stability, reinject the spiked control extract and one or more chromatographic standards at the end of the sample set. Depending on instrument variability and length of sample set, additional spiked control extract or standard injections may be interspersed throughout the sample set.
- vii. Column, Pump, and APCI Interface Care: At the end of set of analyses, flush the column for 30 min with mobile phase D (60/40 acetonitrile/water) at 0.50 mL/min

4. Sample Set
  - a. System Suitability
  - b. External Standard
  - c. Positive Control
  - d. Negative Control
  - e. Intra-laboratory check sample (if needed)
  - f. Samples, up to a maximum of 21
  - g. Reinjection of standard/ positive control

**United States Department of Agriculture  
Food Safety and Inspection Service, Office of Public Health Science**

|                                                                                       |                       |                     |
|---------------------------------------------------------------------------------------|-----------------------|---------------------|
| CLG-FLQ2.01                                                                           |                       | Page 10 of 20       |
| Title: Confirmation of Fluoroquinolone Antibiotics by HPLC Ion Trap Mass Spectrometry |                       |                     |
| Revision: 01                                                                          | Replaces: CLG-FLQ2.00 | Effective: 04/26/21 |

**G. DECISION CRITERIA / CALCULATIONS**

1. Data Processing. Use the QUAL Browser to view total ion current, base ion chromatogram, and/or a total ion chromatogram (TIC) for each drug for each data file. Note retention time of any visible peaks in a drug window. Generate averaged spectra across the retention time window for each drug. This is usually from near the start to near the end of the peak visible in the chromatograms, though a smaller range may be used to avoid a spurious ion spike. Where no peak is visible, use the same settings as in a contemporaneous fortified or positive control extract.
2. Confirmation Criteria:
  - a. Retention times of extract peaks in one or more of the ion chromatograms must match the peak retention time of a contemporaneous (within same analysis set on same day) fortified control extract chromatogram within  $\pm 4\%$ .
  - b. The FLQ peak in the total ion chromatogram (TIC) (see below for ions used for each drug TIC) is present at a S/N ratio of at least 3/1. This is estimated by visual inspection of the TIC.
  - c. The spectrum from the extract must visually match spectra from external standards in the same data set. The base ion must be the same. At least two qualifying ions should be present, readily distinguished from background ions, and have relative abundances comparable to those in the standard. There should be a general absence of nonspecific ions.

Major specific ions for each FLQ are listed below:

| Fluoroquinolone  | Precursor ion(s) | Spectra Range | Base ion | Qualifying Ions |
|------------------|------------------|---------------|----------|-----------------|
| Desethylene cip. | 306              | 80-310        | 306      | 289, 286, 263   |
| Norfloxacin      | 320, 276         | 75-330        | 256      | 257, 233, 219   |
| Ciprofloxacin    | 332, 288         | 75-340        | 268      | 245, 231, 205   |
| Danofloxacin     | 358, 314         | 85-370        | 294      | 283, 245, 219   |
| Desmethyl dano.  | 344              | 90-350        | 344      | 327, 300, 283   |
| Enrofloxacin     | 360, 316         | 85-370        | 245      | 296, 288, 268   |
| Sarafloxacin     | 386, 342         | 90-400        | 322      | 299, 285, 281   |
| Difloxacin       | 400, 356         | 95-410        | 299      | 336, 311, 285,  |

- d. The quality assurance positive and negative control samples confirm and fail to confirm, respectively, for the presence of the appropriate drug.

**United States Department of Agriculture  
Food Safety and Inspection Service, Office of Public Health Science**

|                                                                                       |                       |                     |
|---------------------------------------------------------------------------------------|-----------------------|---------------------|
| CLG-FLQ2.01                                                                           | Page 11 of 20         |                     |
| Title: Confirmation of Fluoroquinolone Antibiotics by HPLC Ion Trap Mass Spectrometry |                       |                     |
| Revision: 01                                                                          | Replaces: CLG-FLQ2.00 | Effective: 04/26/21 |

3. Criteria for Repeating an Analysis:

Note: Sample analyses may be repeated under the following conditions:

- a. The conditions described in G.2.d are not met.
- b. The instrument is suspected to be malfunctioning, as demonstrated by clearly aberrant standard spectra; failure of a calibration check performed shortly after analysis of the sample set; instrumental parameters, especially vacuum readings, outside of normal operating range; or other conditions noted and documented by the analyst.
- c. There is suspected carryover from a previous high concentration sample or standard. In this case, the sample should be reanalyzed after the cause of the carryover has been identified and measures taken to prevent its reoccurrence.
- d. There is strong evidence of FLQ presence, but multiple extraneous ions with relative abundance exceeding that of the FLQ base ion prevent unambiguous confirmation. In this case, it may be appropriate to reanalyze the suspected positive sample together with a chromatographic standard, and negative and positive QA controls.

**H. SAFETY INFORMATION AND PRECAUTIONS**

1. Personal Protective Equipment — Safety glasses, disposable gloves, lab coats.
2. Hazards  
Consult all Safety Data Sheets (SDS) associated with the method
3. Disposal Procedures  
Follow federal, state and local regulations.

**I. QUALITY ASSURANCE PLAN**

1. Performance Standard  
Refer to Section G.2 for Confirmation Criteria.
2. Intralaboratory Check Samples
  - a. Acceptability criteria.  
Refer to I. 1.  
If unacceptable values are obtained, then:
    - i. Investigate following established procedures.
    - ii. Take corrective action as warranted.

**United States Department of Agriculture  
Food Safety and Inspection Service, Office of Public Health Science**

|                                                                                       |                       |                     |
|---------------------------------------------------------------------------------------|-----------------------|---------------------|
| CLG-FLQ2.01                                                                           | Page 12 of 20         |                     |
| Title: Confirmation of Fluoroquinolone Antibiotics by HPLC Ion Trap Mass Spectrometry |                       |                     |
| Revision: 01                                                                          | Replaces: CLG-FLQ2.00 | Effective: 04/26/21 |

**J. APPENDIX**

1. References

Schneider, M. J., Donoghue, D. J. (2002), *J. Chromatogr. B* **780**, 83-92.


2. Chromatograms/spectra

The following chromatograms and spectra are shown on the next 3 pages:

- a. External standards at 25 ppb.
- b. Blank Beef liver.
- c. Beef Liver Recovery at 25 ppb


**United States Department of Agriculture  
Food Safety and Inspection Service, Office of Public Health Science**

|                                                                                       |                       |
|---------------------------------------------------------------------------------------|-----------------------|
| CLG-FLQ2.01                                                                           | Page 13 of 20         |
| Title: Confirmation of Fluoroquinolone Antibiotics by HPLC Ion Trap Mass Spectrometry |                       |
| Revision: 01                                                                          | Replaces: CLG-FLQ2.00 |
|                                                                                       | Effective: 04/26/21   |




**United States Department of Agriculture  
Food Safety and Inspection Service, Office of Public Health Science**

|                                                                                       |                       |                     |
|---------------------------------------------------------------------------------------|-----------------------|---------------------|
| CLG-FLQ2.01                                                                           | Page 14 of 20         |                     |
| Title: Confirmation of Fluoroquinolone Antibiotics by HPLC Ion Trap Mass Spectrometry |                       |                     |
| Revision: 01                                                                          | Replaces: CLG-FLQ2.00 | Effective: 04/26/21 |



**United States Department of Agriculture  
Food Safety and Inspection Service, Office of Public Health Science**









|                                                                                       |                       |
|---------------------------------------------------------------------------------------|-----------------------|
| CLG-FLQ2.01                                                                           | Page 15 of 20         |
| Title: Confirmation of Fluoroquinolone Antibiotics by HPLC Ion Trap Mass Spectrometry |                       |
| Revision: 01                                                                          | Replaces: CLG-FLQ2.00 |
|                                                                                       | Effective: 04/26/21   |



United States Department of Agriculture  
Food Safety and Inspection Service, Office of Public Health Science

|                                                                                       |                       |                     |
|---------------------------------------------------------------------------------------|-----------------------|---------------------|
| CLG-FLQ2.01                                                                           | Page 16 of 20         |                     |
| Title: Confirmation of Fluoroquinolone Antibiotics by HPLC Ion Trap Mass Spectrometry |                       |                     |
| Revision: 01                                                                          | Replaces: CLG-FLQ2.00 | Effective: 04/26/21 |

3. Structures of Fluoroquinolones:






**United States Department of Agriculture  
Food Safety and Inspection Service, Office of Public Health Science**

|                                                                                       |                       |                     |
|---------------------------------------------------------------------------------------|-----------------------|---------------------|
| CLG-FLQ2.01                                                                           | Page 17 of 20         |                     |
| Title: Confirmation of Fluoroquinolone Antibiotics by HPLC Ion Trap Mass Spectrometry |                       |                     |
| Revision: 01                                                                          | Replaces: CLG-FLQ2.00 | Effective: 04/26/21 |

4. Proposed Fragmentation Patterns of Fluoroquinolones:


a. Ciprofloxacin, Sarafloxacin, Norfloxacin, and Enrofloxacin




United States Department of Agriculture  
Food Safety and Inspection Service, Office of Public Health Science

|                                                                                       |                       |
|---------------------------------------------------------------------------------------|-----------------------|
| CLG-FLQ2.01                                                                           | Page 18 of 20         |
| Title: Confirmation of Fluoroquinolone Antibiotics by HPLC Ion Trap Mass Spectrometry |                       |
| Revision: 01                                                                          | Replaces: CLG-FLQ2.00 |
|                                                                                       | Effective: 04/26/21   |

b. Desethylene ciprofloxacin




c. Desmethyl danofloxacin



United States Department of Agriculture  
Food Safety and Inspection Service, Office of Public Health Science

|                                                                                       |                       |                     |
|---------------------------------------------------------------------------------------|-----------------------|---------------------|
| CLG-FLQ2.01                                                                           | Page 19 of 20         |                     |
| Title: Confirmation of Fluoroquinolone Antibiotics by HPLC Ion Trap Mass Spectrometry |                       |                     |
| Revision: 01                                                                          | Replaces: CLG-FLQ2.00 | Effective: 04/26/21 |


d. Danofloxacin



**United States Department of Agriculture  
Food Safety and Inspection Service, Office of Public Health Science**

|                                                                                       |                       |                     |
|---------------------------------------------------------------------------------------|-----------------------|---------------------|
| CLG-FLQ2.01                                                                           | Page 20 of 20         |                     |
| Title: Confirmation of Fluoroquinolone Antibiotics by HPLC Ion Trap Mass Spectrometry |                       |                     |
| Revision: 01                                                                          | Replaces: CLG-FLQ2.00 | Effective: 04/26/21 |

e. Difloxacin



5. Minimum Level of Applicability – 25ppb

**K. APPROVALS AND AUTHORITIES**

1. Approvals on file.
2. Issuing Authority: Director, Laboratory Quality Assurance Staff.