



Food Safety and Inspection Service
U.S. DEPARTMENT OF AGRICULTURE

Listeria monocytogenes Common Findings – Establishment Testing

Brad Webb
Senior Staff Microbiologist
OPPD, RMIS



Inappropriate Sized Sampling Devices


- A 12"x12" area should be sampled, when possible, for FCS and NFCS. If the entire surface area is smaller than 12"x12", then the entire surface should be sampled.
- Smaller sampling devices are not recommended for sampling large areas because they may become saturated.



12 square inches



Smaller devices may be used for small spaces.



Not using a neutralizing broth.

Establishments should utilize a neutralizing broth in FCS and NFCS sampling devices. If not, this could result in false negative results.

USDA FSIS utilizes Dey Engley (DE) Broth because;

- It neutralizes a broad spectrum of antimicrobials, disinfectants, and sanitizers
- DE broth had the greatest efficacy at the time it was adopted.

Other neutralizing broths:

- Letheen Broth
- Neutralizing Buffer
- HiCap Neutralizing Buffer (newest)

Improper storage, shipment, or analyses of samples.

FSIS Guidance;

Samples should be stored and shipped under refrigeration.

- No freezing or excessive heat.
- Refrigeration temps are typically about 2 – 8 C, or 33 – 45 F.



More FSIS Guidance;

- For external labs, limit to 2-3 days from collection to start of analyses.
- For in house lab, initiate testing immediately after collection.
- If not, must support alternative timeframe won't compromise sensitivity.

Sampled Portion vs. Test Portion

- The entire swab or sampling device should be included in the enrichment step.
- Microbiological contamination is non-homogenous.
 - NFCS are typically the original source of contamination of FCS and product.
 - A negative FCS swab only indicates that the specific portion of the FCS that was swabbed was negative at the time it was sampled.
 - A single negative product sample is not proof the entire lot is negative.

Sample Compositing

One sampling device should be utilized per surface location.

If not;

- Potential cross contamination across various surface locations.
- Potential saturation of sampling device.

FSIS recommends that no more than 5 samples be composited;


- Supported by the Lm compliance guide.
- Best if composited by laboratory staff.
- Will lose location specific results.
- Must take corrective actions as if all locations are positive.

Analytical Methods

See [Foodborne Pathogen Test Kits Validated by Independent Organizations \(usda.gov\)](https://www.usda.gov/food-safety/foodborne-pathogen-test-kits-validated-by-independent-organizations).

Establishments and laboratories should choose test kits that are:

1. Validated for testing relevant foods by a:
 - a) Recognized independent body (i.e., AOAC, AFNOR, MicroVal, NordVal),
 - b) U.S. regulatory body (i.e., FSIS MLG or FDA BAM), or
 - c) International Organization for Standardization (ISO) process.



International
Organization for
Standardization

Analytical Methods

- 2) In addition, the validated method should be:
- a) Fit for the intended purpose and application, and
 - b) Performed per the conditions of the validated protocol *by a laboratory that assures the quality of the analytical results.*

Is the lab performing the method as it was validated?

- Look for the matrices that were used in the validation
- Also, look at sample size, enrichment broth(s), times, temperatures, etc.

Submit questions or concerns to askFSIS.

Selection of a Laboratory

Is the lab accredited to ISO 17025 or following FSIS guidance?

Our Selection of a Lab C.G. contains sections on;

- Qualifications of personnel.
- Sample receipt and handling.
- Quality Assurance Management System.
- Method Selection and Implementation.
- Reporting Results.
- Laboratory Assessment Checklist in the Appendix.

Establishment Implementation Issues

- Not listing or testing all FCS sites.
 - Sampling program should include both random and discretionary samples.
 - Should sample the riskiest locations the most often.
- Not testing at the frequency stated.
- Not applying an antimicrobial agent or process AMAP appropriately.
- Not supporting Lm or Sal is not a hazard RLTO for RTE ingredients introduced post-lethality, i.e., spices, sauces, etc.

Implementation Issues

Testing for Lm, but their program states they will test for L. species.


- A product that passes over an FCS which tests positive for Lm is considered adulterated.
- A finding of *Listeria* spp. on an FCS indicates conditions where *Lm* may be present, but the product is not considered adulterated.

Establishments must take corrective actions for all *Listeria* spp. or Lm positives so that product does not become adulterated.



Food Safety and Inspection Service
U.S. DEPARTMENT OF AGRICULTURE

Part 2 – Common Issues with EIAO R_{Lm} or IVT Sampling



EIAO Sampling, Shipping, or Packaging Issues.

- Not check sampling supplies upon receipt, only to realize they're short of supplies when they arrive at the establishment.
- Improper use of gel packs. Frozen gel packs should be removed from the freezer and placed into shipping containers on the day of shipment.
- Breaking establishment procedures due to not reviewing the establishment's operational SSOPs, GMPs, etc., before sampling.

EIAO Sampling, Shipping, or Packaging Issues.

- Allowing an establishment to drastically shorten its RTE lot size to the degree it is no longer representative of routine processing.
- Collecting samples only on the day of packaging, when there was processing, handling, etc., the day before.
- Not understanding the difference between the **sampled lot** vs. a possible **implicated lot** (see page 4 of Directive 10,240.3).

EIAO Sampling, Shipping, or Packaging Issues.

- Examples of issues with aseptic technique;
 - Finger going below the thumb stop on the sponge stick.
 - Sponge contact with the lip of the bag when removing or replacing sponge.
- Collection related mistakes;
 - Not squeezing the broth out of the sponge before removing it from the Whirl-Pak bag.
 - Not squeezing air out of the Whirl-Pak bag before folding down the top.
 - Thinking NFCS samples must be collected near the line.
 - Thinking samples must be collected at pre-op.

EIAO Sampling, Shipping, or Packaging Issues.

Including samples from more than one line or lot within one unit of samples.

- Each unit should be associated with only one line, one production lot, and one 430 Alternative.
- If not, leads to confusion when positives are found.
- Which production lots are implicated by the positive?

Issues with Forms, Seals, Bar-Coded Labels.

- There must be a corresponding form for all samples within each shipping container.
- Each form must be signed.
- Each shipping container should be sealed to ensure sample integrity, as illustrated in Directive 7355.1, revision 3.

DE Broth and Pre-hydrated Sponges

- DE broth can be held at room temperatures for long periods of time, e.g., during sampling.
- Should be refrigerated for long term storage.
- Printed shelf life is dependent on refrigeration.
- Must not be temperature abused.
- The higher the temperature, the faster DE broth degrades.

Questions??

fsis.usda.gov

