



Laboratory Guidebook Notice of Change

Chapter **new**, revised, or archived: MLG 5C Appendix 5.00

Title: PCR Platform Instructions for the real-time PCR detection of Shiga toxin gene and H7 gene in *E. coli* O157:H7

Effective Date: 02/04/2019

Description and purpose of change(s):

This procedure provides instructions for using the ABI® 7500 FAST PCR Platform to detect the *fliCH7* gene encoding the H7 antigen of *E. coli* (including O157:H7) migrating the assay from the SmartCycler® to the ABI® 7500 FAST system. Details include analysis of the data and interpretations for the control results.

**United States Department of Agriculture
Food Safety And Inspection Service, Office of Public Health Science**

MLG 5C Appendix 5.00		Page 1 of 17
Title: PCR Platform Instructions for the real-time PCR detection of Shiga toxin gene and H7 gene in <i>E. coli</i> O157:H7		
Revision: Original	Replaces: NA	Effective: 02/04/2019

Procedure Outline

- A5 5C.1 Introduction
- A5 5C.2 Safety Precautions
- A5 5C.3 Equipment, Reagents, Media and Test Kits
 - A5 5C.3.1 Equipment
 - A5 5C.3.2 Media, Reagents and Cultures
 - A5 5C.3.3 Primers and Probe Sequences
- A5 5C.4 Procedures
 - A5 5C.4.1 Sample preparation
 - A5 5C.4.2 ABI® 7500 FAST PCR Platform Instructions
 - A5 5C.4.3 Preparation of Master Mix
 - A5 5C.4.4 Analyzing Data

A5 5C.1 Introduction

This procedure provides instructions for using the ABI® 7500 FAST PCR Platform to detect the *fliCH7* gene encoding the H7 antigen of *E. coli* (including O157:H7) migrating the assay from the SmartCycler® to the ABI® 7500 FAST system. The real-time PCR assay discerns the shiga toxin (*stx*), intimin (*eae*), *fliCH7* and 16S genes of *E. coli* in a single assay. This assay can be used to confirm these targets in *E. coli* O157:H7 strains. Details include analysis of the data and interpretations for the control results.

A5 5C.2 Safety Precautions

E. coli O157:H7 is categorized as a Biosafety Level 2 human pathogen with a low infectious dose (ingestion of 100 cells can cause disease). CDC guidelines for manipulating Biosafety Level 2 pathogens should be followed whenever live cultures of *E. coli* O157 are used. The use of gloves and eye protection is recommended for all post enrichment viable culture work. Work surfaces should be disinfected prior to and immediately after use. A Class II laminar flow biosafety cabinet is recommended for activities in which infectious aerosols or splashes may be created. The Safety Data Sheets (SDS) must be obtained from the manufacturer for the media, chemicals, reagents and microorganisms used in the analysis. The personnel who will handle the materials should read all SDS.

**United States Department of Agriculture
Food Safety And Inspection Service, Office of Public Health Science**

MLG 5C Appendix 5.00		Page 2 of 17
Title: PCR Platform Instructions for the real-time PCR detection of Shiga toxin gene and H7 gene in <i>E. coli</i> O157:H7		
Revision: Original	Replaces: NA	Effective: 02/04/2019

A5 5C.3 Equipment, Reagents, Media and Test Kits

A5 5C.3.1 Equipment

- a. Incubator ($35 \pm 2^{\circ}\text{C}$)
- b. Heating block ($99 \pm 2^{\circ}\text{C}$) or Thermocycler for DNA template preparation
- c. Mini Vortex mixer
- d. Freezer ($-20 \pm 2^{\circ}\text{C}$)
- e. Clean bench or hood
- f. Applied Biosystems[®] ABI[®] 7500 FAST PCR platform
- g. Applied Biosystems[®] 7500 FAST MicroAMP[™] Fast Optical 96-well Reaction Plate with Barcode (Part# 4346906)
- h. Applied Biosystems[®] 7500 FAST MicroAMP[™] Splash Free 96-Well Base (Part# 4312063)
- i. Applied Biosystems[®] 7500 FAST Optical Adhesive Film (Part# 4311971)

A5 5C.3.2 Media, Reagents and Cultures

- a. Sterile, certified RNase-free, DNase-free molecular grade water (Teknova), or equivalent
- b. 1X Tris-EDTA (TE) Buffer (Teknova Catalog# T0224, or equivalent provider)
- c. Applied Biosystems[®] TaqMan[®] Environmental Master Mix 2.0 (Part# 4396838)
- d. Primers and Probes (IDT Technologies, Coralville, IA)
- e. Tryptic soy agar with 5% sheep blood (SBA) plates, (Remel, Lenexa, KS)
- f. *E. coli* 43894: O157: H7+, stx+, eae+ or equivalent as the positive control

A5 5C.3.3 Primer and Probe Sequences

When receiving lyophilized primers and probes, reconstitute the probe to a stock solution using 1X TE buffer. Working concentrations should be prepared by diluting stocks to the appropriate working concentration. Stocks and working solutions should be stored in the -20°C freezer.

**United States Department of Agriculture
Food Safety And Inspection Service, Office of Public Health Science**

MLG 5C Appendix 5.00		Page 3 of 17
Title: PCR Platform Instructions for the real-time PCR detection of Shiga toxin gene and H7 gene in <i>E. coli</i> O157:H7		
Revision: Original	Replaces: NA	Effective: 02/04/2019

fliC_{H7}

Forward Primer 1068F: 5' TAC CAT CGC AAA AGC AAC TCC 3'
Reverse Primer 1314R: 5' GTC GGC AAC GTT AGT GAT ACC 3'
Probe fliCh7 247: 5'-6-Cy3/ CGG CTG CCG CGA CAT CTT CAA T-IBRQ 3'

stx

Forward Primer: 5'-TTTGTYACTGTSACAGCWGAAGCYTTACS-3'
Reverse Primer: 5'-CCCCAGTTCARWGTRAGRTCMACDTC-3'
Probe stx1: 5'-6-FAM/ZEN/CTGGATGATCTCAGTGGGCGTTCTTATGTAA-IBFQ-3'
Probe stx2: 5'-6-FAM/ZEN/TCGTCAGGCACTGTCTGAAACTGCTCC-IBFQ-3'

eae

Forward Primer: 5'-CATTGATCAGGATTTTTCTGGTGATA-3'
Reverse Primer: 5'-CTCATGCGGAAATAGCCGTTM-3'
Probe: 5'-MAX/ATAGTCTCGCCAGTATTCGCCACCAATACC-IBFQ-3'

16S rRNA Internal Control

Forward Primer: 5' CCT CTT GCC ATC GGA TGT G-3'
Reverse Primer: 5' GGC TGG TCA TCC TCT CAG ACC-3'
Probe: 5' TYE665/GTG GGG TAA CGG CTC ACC TAG GCG AC-IBRQ-3'

A5 5C.4 Procedures

A5 5C.4.1 Sample preparation

Note: A template from a 16-24 h isolate on SBA should be used for the PCR assay but the template can be prepared from a selective medium, if desired.

- a. For each sample and the control (H7 and *stx* positive *E. coli*), streak a SBA plate to obtain isolated colonies and incubate at 35 ± 2°C for 16-24 hrs.
- b. For each sample and control, dispense 50 ± 1µl of sterile molecular grade water into a labeled sterile micro centrifuge tube (for use on heating block) or thin-walled PCR tube (for use on thermocycler).
- c. Transfer one isolated colony from Step 1 into each labeled tube. One tube will have no colony (only water) as a blank.
- d. Heat the tubes in a thermocycler or heating block for 10 minutes at 99 ± 2°C. These samples will be templates for PCR.

**United States Department of Agriculture
Food Safety And Inspection Service, Office of Public Health Science**

MLG 5C Appendix 5.00	Page 4 of 17	
Title: PCR Platform Instructions for the real-time PCR detection of Shiga toxin gene and H7 gene in <i>E. coli</i> O157:H7		
Revision: Original	Replaces: NA	Effective: 02/04/2019

- e. Keep tubes at 2-8°C until ready to use.

A5 5C.4.2 ABI® 7500 FAST PCR Platform Instructions

- a. Turn on the ABI® 7500 FAST machine and computer.
b. If you are just looking at the data without the machine on, select “Continue without connecting to machine” to proceed.
c. Open the ABI® software and log on as guest.
d. Choose the advanced set-up button on the left.



- e. Enter the name of the file according to the laboratory-naming scheme (e.g. MF12345).
f. Ensure that items appear as below [instrument: 7500 FAST; experiment: Quantitation-Standard Curve, reagents-Taqman Reagents), except change the ramp speed from FAST to Standard chemistry.

Which instrument are you using to run the experiment?

7500 (96 Wells) 7500 Fast (96 Wells)

Set up, run, and analyze an experiment using a fast cycling 5-color, 96-well system.

What type of experiment do you want to set up?

Quantitation - Standard Curve Quantitation - Relative Standard Curve Quantitation - Comparative Ct ($\Delta\Delta CT$)

Melt Curve Genotyping Presence/Absence

Use standards to determine the absolute quantity of target nucleic acid sequence in samples.

Which reagents do you want to use to detect the target sequence?

TaqMan® Reagents SYBR® Green Reagents Other

The PCR reactions contain primers designed to amplify the target sequence and a TaqMan® probe designed to detect amplification of the target sequence.

Which ramp speed do you want to use in the instrument run?

Standard (~ 2 hours to complete a run) Fast (~ 40 minutes to complete a run)

For optimal results with the standard ramp speed, Applied Biosystems recommends using standard reagents for your PCR reactions.

- g. Click on the plate set-up tab on the left



**United States Department of Agriculture
Food Safety And Inspection Service, Office of Public Health Science**

MLG 5C Appendix 5.00		Page 5 of 17
Title: PCR Platform Instructions for the real-time PCR detection of Shiga toxin gene and H7 gene in <i>E. coli</i> O157:H7		
Revision: Original	Replaces: NA	Effective: 02/04/2019

h. Click on “Define Targets and Samples” tab. Enter the targets for the assays listed in Table 1.





Table 1. *Salmonella* confirmatory PCR targets

Target	Dye	Quencher
fliCH7	Cy3	None
stx	FAM	None
Eae	Joe	None
16S	Cy5	None

Under “Define Samples”, click on the number of samples button and add samples until you have enough for all your assays, including control samples (positive control and NTC).



**United States Department of Agriculture
Food Safety And Inspection Service, Office of Public Health Science**

MLG 5C Appendix 5.00		Page 6 of 17
Title: PCR Platform Instructions for the real-time PCR detection of Shiga toxin gene and H7 gene in <i>E. coli</i> O157:H7		
Revision: Original	Replaces: NA	Effective: 02/04/2019

- i. After you change the sample name (or do this after the run) and choose any color changes (optional) you want, click on the “Assign Target and Samples” tab.



Highlight the well that you want to place the sample.

The screenshot shows the 'Assign Targets and Samples' software interface. It has two tabs: 'Define Targets and Samples' and 'Assign Targets and Samples'. The 'Assign Targets and Samples' tab is active. On the left, there are instructions and two main sections: 'Assign target(s) to the selected wells.' and 'Assign sample(s) to the selected wells.'. The 'Assign target(s) to the selected wells.' section has a table with columns 'Assign', 'Target', 'Task', and 'Quant'. The 'Assign' column has checkboxes. The 'Target' column lists 'fliH7', 'stx', 'eae', and '16S'. The 'Task' column has buttons for 'U' (Unknown), 'S' (Standard), and 'N' (Negative Control). The 'Assign sample(s) to the selected wells.' section has a table with columns 'Assign' and 'Sample'. The 'Assign' column has checkboxes. The 'Sample' column lists 'Sample 1', 'Sample 2', and 'Sample 3'. On the right, there is a 'View Plate Layout' section with a grid of wells. The grid has columns 1-7 and rows A-F. The well at row A, column 1 is highlighted with a red circle. There are also buttons for 'Show in Wells' and 'View Legend'.

Select the sample from the assign samples box by checking the box. Repeat until all samples are placed in the correct locations.

**United States Department of Agriculture
Food Safety And Inspection Service, Office of Public Health Science**

MLG 5C Appendix 5.00		Page 7 of 17
Title: PCR Platform Instructions for the real-time PCR detection of Shiga toxin gene and H7 gene in <i>E. coli</i> O157:H7		
Revision: Original	Replaces: NA	Effective: 02/04/2019

Define Targets and Samples
Assign Targets and Samples

I Instructions: To set up standards: Click "Define and Set Up Standards."
 To set up unknowns: Select wells, assign target(s), select "U" (Unknown) as the task for each target.
 To set up negative controls: Select wells, assign target(s), then select "N" (Negative Control) as the task for each target.

Assign target(s) to the selected wells.

Assign	Target	Task	Quant
<input type="checkbox"/>	fliH7	<input type="checkbox"/> U <input type="checkbox"/> S <input type="checkbox"/> N	
<input type="checkbox"/>	stx	<input type="checkbox"/> U <input type="checkbox"/> S <input type="checkbox"/> N	
<input type="checkbox"/>	eae	<input type="checkbox"/> U <input type="checkbox"/> S <input type="checkbox"/> N	
<input type="checkbox"/>	16S	<input type="checkbox"/> U <input type="checkbox"/> S <input type="checkbox"/> N	

Mixed
 U Unknown
 S Standard
 N Negative Control

Define and Set Up Standards

View Plate Layout View Well

Select Wells View

Show in Wells ▼ View Legend

	1	2	3	4
A	Sample 1			
B				
C				
D				
E				
F				
G				
H				

Assign sample(s) to the selected wells.

Assign	Sample
<input checked="" type="checkbox"/>	Sample 1
<input type="checkbox"/>	Sample 2
<input type="checkbox"/>	Sample 3

Assign sample(s) of selected well(s) to biological group(s).

Assign	Biological Group
<input type="checkbox"/>	

Select the dye to use as the passive reference.

ROX ▼

**United States Department of Agriculture
Food Safety And Inspection Service, Office of Public Health Science**

MLG 5C Appendix 5.00		Page 8 of 17
Title: PCR Platform Instructions for the real-time PCR detection of Shiga toxin gene and H7 gene in <i>E. coli</i> O157:H7		
Revision: Original	Replaces: NA	Effective: 02/04/2019

Highlight all wells that contain samples, hold down the CTRL button to add multiple wells.

Define Targets and Samples
Assign Targets and Samples

Instructions:
 To set up standards: Click "Define and Set Up Standards."
 To set up unknowns: Select wells, assign target(s), select "U" (Unknown) as the task for each target assign
 To set up negative controls: Select wells, assign target(s), then select "N" (Negative Control) as the task for

Assign target(s) to the selected wells.

Assign	Target	Task	Quant
<input type="checkbox"/>	fliH7	<input type="checkbox"/> U <input type="checkbox"/> S <input type="checkbox"/> N	
<input type="checkbox"/>	stx	<input type="checkbox"/> U <input type="checkbox"/> S <input type="checkbox"/> N	
<input type="checkbox"/>	eae	<input type="checkbox"/> U <input type="checkbox"/> S <input type="checkbox"/> N	
<input type="checkbox"/>	16S	<input type="checkbox"/> U <input type="checkbox"/> S <input type="checkbox"/> N	

* Mixed U Unknown S Standard N Negative Contr

Define and Set Up Standards

Assign sample(s) to the selected wells.

Assign	Sample
<input type="checkbox"/>	Sample 1
<input type="checkbox"/>	Sample 2
<input type="checkbox"/>	Sample 3

Assign sample(s) of selected well(s) to biological

Assign	Biological Group
<input type="checkbox"/>	

Select the dye to use as the passive reference.

ROX

View Plate Layout

View Well Table

Select Wells With:

	1	2	3	4	5	6
A	Sample 1					
B	Sample 2					
C	Sample 3					
D	Sample 4					
E	Sample 5					
F	Sample 6					
G	Sample 7					
H	Sample 8					

After you assign all of your samples to the wells, select the targets (fliH7, stx, eae and 16S) you want from the assign targets box. It will bring up 4 colors in each well box which is why it is best to put sample names before this step so that you can see them.

**United States Department of Agriculture
Food Safety And Inspection Service, Office of Public Health Science**

MLG 5C Appendix 5.00		Page 9 of 17
Title: PCR Platform Instructions for the real-time PCR detection of Shiga toxin gene and H7 gene in <i>E. coli</i> O157:H7		
Revision: Original	Replaces: NA	Effective: 02/04/2019

Define Targets and Samples
Assign Targets and Samples

I Instructions: To set up standards: Click "Define and Set Up Standards."
 To set up unknowns: Select wells, assign target(s), select "U" (Unknown) as
 To set up negative controls: Select wells, assign target(s), then select "N" (N

Assign target(s) to the selected wells.

Assign	Target	Task	Quant
<input checked="" type="checkbox"/>	fliH7	U S N	
<input checked="" type="checkbox"/>	stx	U S N	
<input checked="" type="checkbox"/>	eae	U S N	
<input checked="" type="checkbox"/>	16S	U S N	

* Mixed U Unknown S Standard N Negative Con

Define and Set Up Standards

Assign sample(s) to the selected wells.

Assign	Sample
<input type="checkbox"/>	Sample 1
<input type="checkbox"/>	Sample 2
<input type="checkbox"/>	Sample 3

Assign sample(s) of selected well(s) to biological

Assign	Biological Group
<input type="checkbox"/>	

Select the dye to use as the passive reference.

ROX

View Plate Layout

Show in Wells ▼

	1	2
A	U eae U fliH7	
B	U eae U fliH7	
C	U eae U fliH7	
D	U eae U fliH7	
E	U eae U fliH7	
F	U eae U fliH7	
G	U eae U fliH7	
H	U eae U fliH7	

**United States Department of Agriculture
Food Safety And Inspection Service, Office of Public Health Science**

MLG 5C Appendix 5.00		Page 10 of 17
Title: PCR Platform Instructions for the real-time PCR detection of Shiga toxin gene and H7 gene in <i>E. coli</i> O157:H7		
Revision: Original	Replaces: NA	Effective: 02/04/2019

In order to see all of the targets and sample names in each well you may stretch the wells using the arrows. This will allow all the contents of the well to be seen.

Define Targets and Samples | **Assign Targets and Samples**

Instructions:
 To set up standards: Click "Define and Set Up Standards."
 To set up unknowns: Select wells, assign target(s), select "U" (Unknown) as th
 To set up negative controls: Select wells, assign target(s), then select "N" (Neg

Assign target(s) to the selected wells.

Assign	Target	Task	Quant
<input checked="" type="checkbox"/>	fliH7	<input type="checkbox"/> U <input type="checkbox"/> S <input type="checkbox"/> N	
<input checked="" type="checkbox"/>	stx	<input type="checkbox"/> U <input type="checkbox"/> S <input type="checkbox"/> N	
<input checked="" type="checkbox"/>	eae	<input type="checkbox"/> U <input type="checkbox"/> S <input type="checkbox"/> N	
<input checked="" type="checkbox"/>	16S	<input type="checkbox"/> U <input type="checkbox"/> S <input type="checkbox"/> N	

Mixed Unknown Standard Negative Con

Define and Set Up Standards

Assign sample(s) to the selected wells.

Assign	Sample
<input type="checkbox"/>	Sample 1
<input type="checkbox"/>	Sample 2
<input type="checkbox"/>	Sample 3

Assign sample(s) of selected well(s) to biological

Assign	Biological Group
<input type="checkbox"/>	

Select the dye to use as the passive reference.

ROX

View Plate Layout

Show in Wells

1 2

A
Sample 1
U 16S
U eae
U fliH7
U stx

B
Sample 2
U 16S
U eae
U fliH7
U stx

C
Sample 3
U 16S
U eae
U fliH7
U stx


D
Sample 4
U 16S
U eae
U fliH7
U stx

Wells: U 8 Unknown S 0


- j. Make sure the little box that is at the very bottom of the page says passive reference dye set to "ROX".

**United States Department of Agriculture
Food Safety And Inspection Service, Office of Public Health Science**

MLG 5C Appendix 5.00		Page 11 of 17
Title: PCR Platform Instructions for the real-time PCR detection of Shiga toxin gene and H7 gene in <i>E. coli</i> O157:H7		
Revision: Original	Replaces: NA	Effective: 02/04/2019



- k. **(OPTIONAL)** If running a standard curve, click the orange S box before assigning the wells. Otherwise, all samples are unknown, which is the default so that parameter should not require checking.



- l. Click on the left tab to set-up run method.



Remove the first holding step by right clicking that box and deleting.
Change the temperatures and cycles to the parameters defined in Table 2.


Table 2. PCR assay parameters

Number of Cycles	Temperature	Time
1 (Holding Stage)	95°C	10 minutes
35 (Cycling Stages)	95°C	15 seconds
	59°C	60 seconds

**United States Department of Agriculture
Food Safety And Inspection Service, Office of Public Health Science**

MLG 5C Appendix 5.00		Page 12 of 17
Title: PCR Platform Instructions for the real-time PCR detection of Shiga toxin gene and H7 gene in <i>E. coli</i> O157:H7		
Revision: Original	Replaces: NA	Effective: 02/04/2019

To change parameters, click on the temperature and change to the appropriate setting for that step. Next, click on the time and change to the appropriate setting for that step.



- m. At the top of the graph, change reaction volume per well to 25 µL.
Reaction Volume Per Well µL
- n. Prepare your PCR assay plate at this step.

A5 5C.4.3 Preparation of Master Mix

Perform this procedure in a designated pre-PCR area (clean bench or hood). Always wear gloves while working in this area. Keep the pipettors in this area only for PCR preparation. Keep all reagents and Master Mix at 2-8°C at all times.

**United States Department of Agriculture
Food Safety And Inspection Service, Office of Public Health Science**

MLG 5C Appendix 5.00		Page 13 of 17
Title: PCR Platform Instructions for the real-time PCR detection of Shiga toxin gene and H7 gene in <i>E. coli</i> O157:H7		
Revision: Original	Replaces: NA	Effective: 02/04/2019

Probes are light sensitive, so it is necessary to limit exposure to light as much as possible.

Use two controls for each assay:

- *E. coli* O157:H7, toxin positive
 - No Template Control (NTC)– Molecular Grade water
- a. Calculate the number of samples plus controls to be tested (plus at least one to compensate for pipettor error). Calculate the volume of each primer, probe, ABI® Environmental Mastermix and water to be added to make the master mix by following the chart below.

Reaction vol (µL)	25
Number of rxns	1

Component	µL
ABI® Mastermix	12.50 µL x # of reactions
primer 16S F (20uM working)	0.20 µL x # of reactions
primer 16S R (20uM working)	0.20 µL x # of reactions
fliCH7 F (100uM working)	0.25 µL x # of reactions
fliCH7 R (100uM working)	0.25 µL x # of reactions
Stx F (100uM working)	0.31 µL x # of reactions
Stx R (100uM working)	0.31 µL x # of reactions
Eae F (100um working)	0.25 µL x # of reactions
Eae R (100um working)	0.25 µL x # of reactions

**United States Department of Agriculture
Food Safety And Inspection Service, Office of Public Health Science**

MLG 5C Appendix 5.00		Page 14 of 17
Title: PCR Platform Instructions for the real-time PCR detection of Shiga toxin gene and H7 gene in <i>E. coli</i> O157:H7		
Revision: Original	Replaces: NA	Effective: 02/04/2019

16S P (5uM working)	0.50 µL x # of reactions
fliCH7 P (5uM working)	1.0 µL x # of reactions
Stx1 P (5uM working)	1.25 µL x # of reactions
Stx2 P (5uM working)	1.25 µL x # of reactions
Eae P (5uM working)	1.0 µL x # of reactions
dH2O	0.475 µL x # of reactions
Template DNA	5uL

- b. Add each reagent accordingly to a labeled microcentrifuge tube.
- c. Mix by inversion or vortexing. Briefly centrifuge the mixture to remove any liquid from the lid.
- d. Aliquot 20.0 µl of the master mix into the designated tubes or wells of the 96 well plate.
- e. Add 5.0 µl of template to the appropriate tubes or wells.
- f. Close the tubes.
- g. Centrifuge the tubes or 96 well plate using a pulse spin setting. Ensure that there are no bubbles present at the bottom of the wells.
- h. After you have the mastermix and template added to the appropriate wells, seal the plate with the MicroAmp Optical Adhesive Film and centrifuge the plate briefly to ensure no bubbles are at the bottom of each well.
- i. On the ABI® 7500 FAST instrument, open the holder for the plate by pressing the indentation next to the power button. Change the holder from strips to 96-well plate if needed. If you use the strips, you MUST put extra strips in the holder on the opposite side of the tube holder and they must be pushed into the grooves perfectly to avoid crushing your tubes. If you are using more than three columns of reactions (24 reactions), it is more economical to use the 96-well plates and seals.
- j. Set the plate or strips into the holder and re-close using the indentation on the right.
- k. In the software interface, the temperature cycling graph should still be up. Click the large Green button labeled “Start Run”. Store the experiment run


**United States Department of Agriculture
Food Safety And Inspection Service, Office of Public Health Science**

MLG 5C Appendix 5.00		Page 15 of 17
Title: PCR Platform Instructions for the real-time PCR detection of Shiga toxin gene and H7 gene in <i>E. coli</i> O157:H7		
Revision: Original	Replaces: NA	Effective: 02/04/2019

data in a location accessible to the laboratory. If you do not want to make any changes to where the file is stored, click open and the machine will start running.




- l. The start run button will turn from green to red and it will show the amount of time left and what step is running in the program.
- m. Do not edit or click anything while it is running. There should be a graph showing all the raw data as the program proceeds.
- n. After it is finished, take out your plate, close the door and turn off the machine.
- o. Press “save” at the top of the screen. If you want to analyze the data at a location with the ABI® 7500 Software, move the .eds file to the designated folder on the server and open within the program. Otherwise, click the green analyze button at the top.



A5 5C.4.4 Analyzing Data

- a. To see the amplification curves, choose target: invA and/or 16S at the bottom left options panel. It starts by showing all three in the log phase. Change to linear phase.




- b. To see the threshold automatically determined by the software, check the show threshold box. Auto will be checked since the threshold has been set automatically. A thick line with the threshold for the selected target should come up. To see the threshold line and value, leave as is with all three and it will show three thresholds or click on a single target using the drop-down box.
- c. To show curves of the selected samples, click on a well in the 96-well plate diagram. Use the ctrl button to click multiple wells simultaneously. **NOTE:**

**United States Department of Agriculture
Food Safety And Inspection Service, Office of Public Health Science**

MLG 5C Appendix 5.00		Page 16 of 17
Title: PCR Platform Instructions for the real-time PCR detection of Shiga toxin gene and H7 gene in <i>E. coli</i> O157:H7		
Revision: Original	Replaces: NA	Effective: 02/04/2019

Any time you make a change (e.g. sample name or location), click the green analyze button for the changes to take effect.



- d. To view the C_T values, click from the “View Plate Layout” tab to the “View Well Table” tab. If you do not see any C_T values, click the analyze button again.
- e. Choose “Group by Target Name” to view results by gene target.




- f. To **export** results (make sure you pressed “Analyze” button before this step), click “Export” on the top bar. To export the data for all samples, **make sure that all samples are selected on the graph**. Choose the destination for the file and click “Customize Export” to check which boxes you want exported. To simplify data analyses, select the following items from “Select Results Content”: Well, Sample Name, Target Name, Reporter, and C_T.

Well	Sample ...	Reporter	Ct	Target N...
A1	1	CY5	14.654549	16S
A1	1	FAM	15.953678	invA
B1	2	CY5	14.952367	16S
B1	2	FAM	15.755151	invA
C1	3	CY5	14.836021	16S
C1	3	FAM	15.257235	invA
D1	4	CY5	24.409437	16S
D1	4	FAM	24.3158	invA
E1	5	CY5	15.096321	16S
E1	5	FAM	15.643331	invA
F1	6	CY5	14.138546	16S
F1	6	FAM	13.537354	invA
G1	7	CY5	12.706549	16S
G1	7	FAM	Undetermi...	invA
H1	pos	CY5	17.212683	16S
H1	pos	FAM	17.10175	invA

**United States Department of Agriculture
Food Safety And Inspection Service, Office of Public Health Science**

MLG 5C Appendix 5.00		Page 17 of 17
Title: PCR Platform Instructions for the real-time PCR detection of Shiga toxin gene and H7 gene in <i>E. coli</i> O157:H7		
Revision: Original	Replaces: NA	Effective: 02/04/2019

- g. **(OPTIONAL)** If you want to save the **picture** of the graph click the button to the right of the print button at the top of the graph and save as a jpg.



- h. Data can be analyzed directly from the ABI[®] 7500 FAST machine or from the exported file. A reported C_T value (e.g. 24.5 for *fliH7*) for a target indicates that the sample is positive for that target.