

**United States Department of Agriculture
Food Safety and Inspection Service, Office of Public Health Science**

CLG-AVR.02		Page 1 of 17
Title: Determination of Ivermectin, Doramectin, and Moxidectin by HPLC		
Revision: 02	Replaces: CLG-AVR.01	Effective: 12/22/2006

Contents

A.	INTRODUCTION	2
B.	EQUIPMENT	5
C.	REAGENTS AND SOLUTIONS	6
D.	STANDARDS	7
E.	SAMPLE PREPARATION	8
F.	ANALYTICAL PROCEDURE.....	8
G.	CALCULATIONS	11
H.	SAFETY INFORMATION AND PRECAUTIONS.....	12
I.	QUALITY ASSURANCE PLAN	13
J.	WORKSHEET	15
K.	APPENDIX	17
L.	APPROVALS AND AUTHORITIES	17

**United States Department of Agriculture
Food Safety and Inspection Service, Office of Public Health Science**

CLG-AVR.02		Page 2 of 17
Title: Determination of Ivermectin, Doramectin, and Moxidectin by HPLC		
Revision: 02	Replaces: CLG-AVR.01	Effective: 12/22/2006

A. INTRODUCTION

1. Theory

Moxidectin, ivermectin and doramectin are potent anthelmintics used in food animals to control parasitic infections. Moxidectin, ivermectin and doramectin are extracted from tissue with acetonitrile; extraneous substances are removed using alumina chromatographic cleanup. The analytes are determined by HPLC after formation of fluorescent derivatives with trifluoroacetic anhydride/1-methylimidazole.


2. Applicability

The method is applicable to determination of Moxidectin, Ivermectin, and Doramectin in liver and muscle of bovine, ovine, porcine, caprine and equine species at ≥ 7.5 ppb.


**United States Department of Agriculture
Food Safety and Inspection Service, Office of Public Health Science**

CLG-AVR.02	Page 3 of 17
Title: Determination of Ivermectin, Doramectin, and Moxidectin by HPLC	
Revision: 02	Replaces: CLG-AVR.01 Effective: 12/22/2006


3. Structures



Doramectin




Doramectin Derivative



Ivermectin

R = C₂H₅ for H₂B_{1A}
R = CH₃ for H₂B_{1B}




Ivermectin Derivative

R = C₂H₅ for H₂B_{1A}
R = CH₃ for H₂B_{1B}

**United States Department of Agriculture
Food Safety and Inspection Service, Office of Public Health Science**


CLG-AVR.02		Page 4 of 17
Title: Determination of Ivermectin, Doramectin, and Moxidectin by HPLC		
Revision: 02	Replaces: CLG-AVR.01	Effective: 12/22/2006

3. Structures (cont.)



Abamectin

R = C₂H₅ for H₂B_{1A}
R = CH₃ for H₂B_{1B}




Abamectin Derivative

R = C₂H₅ for H₂B_{1A}
R = CH₃ for H₂B_{1B}

**United States Department of Agriculture
Food Safety and Inspection Service, Office of Public Health Science**

CLG-AVR.02		Page 5 of 17
Title: Determination of Ivermectin, Doramectin, and Moxidectin by HPLC		
Revision: 02	Replaces: CLG-AVR.01	Effective: 12/22/2006

3. Structures (cont.)



Moxidectin

B. EQUIPMENT

Note: Equivalent apparatus and/or instrumentation may be substituted.

1. Apparatus

- a. N-EVAP - Model 112, Organomation Assoc. Inc.
- b. Centrifuge - Sorvall model T-6000B, DuPont Co.
- c. Mechanical shaker - Eberbach model 610 equipped with shaker box model 6040, Thomas Scientific.
- d. Vortex mixer - Fisher Scientific.
- e. Extraction columns - Fisher Scientific Prep Sep-R (empty), Cat. No. P449R, Fisher Scientific.
- f. 50 mL screw cap centrifuge tubes - Cat. No. 05-558-12B, Fisher Scientific.
- g. 50 mL polypropylene centrifuge tubes - Cat. No. 222-3937-G80, Evergreen Scientific International Inc.
- h. EDP Plus Micropipet - Rainin Instruments Inc.
- i. Eppendorf Pipettor - Cat. No. 4789, Brinkman Instrument Inc.
- j. Eppendorf Combitips - 5.0 mL and 12.5 mL, Brinkman Instrument Inc.
- k. SPE Cartridges - Place a small silanized glass wool plug into the neck of a 5.75" disposable transfer pipet. Add 0.1 ± 0.01 g C18 bulk packing material into the

**United States Department of Agriculture
Food Safety and Inspection Service, Office of Public Health Science**

CLG-AVR.02		Page 6 of 17
Title: Determination of Ivermectin, Doramectin, and Moxidectin by HPLC		
Revision: 02	Replaces: CLG-AVR.01	Effective: 12/22/2006

disposable pipet. Tap gently to settle.

- I. Glass test tubes - 16 x 100 mm, 20 x 150 mm, and 12 x 75 mm, Fisher Scientific.

2. Instrumentation

a. Liquid Chromatography System

- i. Agilent 1100 Series Quaternary Pump, Agilent Technologies, Inc.
- ii. Agilent 1100 Series Autosampler, Agilent Technologies, Inc.
- iii. Agilent 1100 Series Vacuum Degasser, Agilent Technologies, Inc.
- iv. Agilent 1100 Series Fluorescence detector, Agilent Technologies, Inc.
- v. Agilent 1100 Series Column Compartment, Agilent Technologies, Inc.
- vi. Zorbax ODS 4.6 mm x 15 cm C18 analytical column, Agilent Technologies, Inc.
- vii. Brownlee Columns Spheri-5 RP-18, 30 mm x 4.6 mm guard column, 5 micron particle size, Perkin Elmer.

C. REAGENTS AND SOLUTIONS

Note: Equivalent reagents and solutions may be substituted if necessary.

1. Reagents

- a. Acetonitrile LC Grade
- b. Alumina-Neutral type WN-3, Activity grade 1, Sigma Chemical Co. Dry at 135 ± 5 °C for at least 24 hours prior to use.

Prepare deactivated alumina for column chromatography. Alumina should be 12% deactivated. For example: Add 24 g distilled/deionized water to 176 g alumina. Mix by shaking until there are no visible lumps. Store deactivated alumina at room temperature in a tightly closed container. Use within one week after opening.

Prepare alumina columns by weighing 2.0 ± 0.2 grams of deactivated alumina into an empty Prep-Sep column.
- c. 1-Methylimidazole - redistilled (99+ %), Cat. No. 33,609-2, Aldrich Chemical Co.
- d. Trifluoroacetic anhydride (TFAA) - (99+ %), Cat. No. 10,623-2, Aldrich Chemical Co.
- e. Methanol - LC Grade.
- f. Water - HPLC grade.

**United States Department of Agriculture
Food Safety and Inspection Service, Office of Public Health Science**

CLG-AVR.02		Page 7 of 17
Title: Determination of Ivermectin, Doramectin, and Moxidectin by HPLC		
Revision: 02	Replaces: CLG-AVR.01	Effective: 12/22/2006

2. Solutions

- a. 1-Methylimidazole - 1:1 v/v 1-methylimidazole/acetonitrile:
Add 1 part acetonitrile to 1 part 1-methylimidazole.
For forty samples, add 5 mL of acetonitrile to 5 ml of 1-methylimidazole and mix.
- b. TFAA - 1:1 v/v trifluoroacetic anhydride/acetonitrile:
Add 1 part of acetonitrile to 1 part trifluoroacetic anhydride.
For forty samples, add 5 mL of acetonitrile and 5 mL of trifluoroacetic anhydride and mix.

D. STANDARDS

1. Source

- a. Ivermectin standard, Cat. No. L-640,471-076P004, Merck Manufacturing Division, West Point, PA.
- b. Abamectin standard, Cat. No. L-676,863-038A003, Merck Manufacturing Division, West Point, PA.
- c. Doramectin standard, Cat. No. UK-67,994, Pfizer, Groton, CT.
- d. Moxidectin standard, Cat. No. 301423, Fort Dodge Animal Health.

2. Preparation

- a. Stock solution:
Follow manufacturer's instructions accompanying standards to obtain a stock solution of approximately 125 ± 1 µg/mL in acetonitrile.
- b. Fortification solution:
Make a 1:250 dilution to obtain a fortification solution of 0.5 µg/mL in acetonitrile.
Add 1 mL each of doramectin, ivermectin and moxidectin stock solutions to a 250 mL flask and bring to volume with acetonitrile and mix.

3. Storage and Stability (if not included with preparation)

- a. Store stock solution in freezer < -10 °C.
- b. Fortification solutions may be stored at room temperature.
- c. Stability:

Fortification solution: 90 days.

Stock solution: 1 year.

**United States Department of Agriculture
Food Safety and Inspection Service, Office of Public Health Science**

CLG-AVR.02		Page 8 of 17
Title: Determination of Ivermectin, Doramectin, and Moxidectin by HPLC		
Revision: 02	Replaces: CLG-AVR.01	Effective: 12/22/2006

E. SAMPLE PREPARATION

Process samples of liver and muscle until homogeneous. All samples are stored refrigerated or frozen until analyzed.

Note: Remove excess fat from tissue before homogenization.

F. ANALYTICAL PROCEDURE

1. Sample extraction

- a. Weigh 2.5 ± 0.2 g ground tissue into a 50 mL polypropylene centrifuge tube.
- b. Add 8 mL acetonitrile and vortex for 30 sec.
- c. Fortify each sample with 150 μ L of 0.5 μ g/mL (equivalent to 30 ppb) abamectin internal standard solution. Fortify moxidectin, ivermectin and doramectin recoveries with appropriate volume of 0.5 μ g/mL fortification solution depending on species analyzed.
- d. Centrifuge for 3 min at 1500 RPM.
- e. Pour acetonitrile eluent through deactivated alumina column and collect eluate in a 50 mL glass centrifuge tube.
- f. Repeat extraction with additional 8 mL acetonitrile, centrifuge and decant through alumina column combining eluents.
- g. Evaporate acetonitrile under a gentle stream of dry nitrogen or dry air at 65 °C.
- h. Reconstitute the dried sample from step (g) using 0.5 mL acetonitrile. Vortex to mix.
- i. For muscle tissue - Add 2 mL acetonitrile and proceed to step (m).

For liver tissue - Prepare SPE C18 cartridge by placing a small silanized glass wool plug into the neck of a 5.75" disposable pipet. Add 0.1 g of C18 bulk packing material into the disposable pipet on top of the glass wool. Tap gently to settle C18. This SPE clean-up step eliminates co-extractants that interfere with quantitation.

- j. Pre-wet the SPE cartridge with 1.0 mL acetonitrile. Discard the wash.
- k. Load the 0.5 mL sample from step (h) onto the wet SPE cartridge. **(The columns must not dry out at any time or the Avermectin recoveries will be low).** Collect the eluent in a test tube.
- l. Add 2 mL of acetonitrile to the sample tube and mix. Add to the SPE column. Collect the eluent in the same container as the initial 0.5 mL eluent.

Note: Properly dispose of the SPE column.

- m. Add 200 μ L 1-methylimidazole/acetonitrile reagent to the eluent and vortex at

**United States Department of Agriculture
Food Safety and Inspection Service, Office of Public Health Science**

CLG-AVR.02		Page 9 of 17
Title: Determination of Ivermectin, Doramectin, and Moxidectin by HPLC		
Revision: 02	Replaces: CLG-AVR.01	Effective: 12/22/2006

least 10 sec.

- n. Add 200 µL TFAA/acetonitrile reagent and vortex at least 10 sec.
- o. Allow the sample to derivatize in the dark for a minimum of 15 minutes before HPLC analysis.

Note: Derivatized samples decompose on exposure to strong light.

- p. Inject on HPLC.

2. Instrumental Settings


Note: System may be adjusted to insure optimum response.

- a. Mobile phase: 3:97 v/v water/methanol.
- b. Flow rate: 1.8 mL /min.
- c. Column temp.: 30 °C.
- d. Injection vol.: 50 µL - As determined by detector/integrator conditions.
- e. Run time: 15 min.
- f. Detector settings:
 - i. Excitation wavelength 365 ± 20 nm.
 - ii. Emission wavelength 465 ± 20 nm.


**United States Department of Agriculture
Food Safety and Inspection Service, Office of Public Health Science**

CLG-AVR.02		Page 10 of 17
Title: Determination of Ivermectin, Doramectin, and Moxidectin by HPLC		
Revision: 02	Replaces: CLG-AVR.01	Effective: 12/22/2006

3. Sample Chromatograms



a. External Standard of Moxidectin (7.5 ppb), Abamectin (15 ppb), Doramectin (7.5 ppb), and Ivermectin (7.5 ppb).



b. Beef Liver fortified with Moxidectin (7.5 ppb), Abamectin (15 ppb), Doramectin (7.5 ppb), and Ivermectin (7.5 ppb).

**United States Department of Agriculture
Food Safety and Inspection Service, Office of Public Health Science**

CLG-AVR.02		Page 11 of 17
Title: Determination of Ivermectin, Doramectin, and Moxidectin by HPLC		
Revision: 02	Replaces: CLG-AVR.01	Effective: 12/22/2006

G. CALCULATIONS

Quantitation is performed by measuring peak area. Each set is accompanied by an external standard curve at 0, 7.5, 15, 30, and 60 ppb.

Measure peak area of abamectin, ivermectin, moxidectin, and doramectin peaks in the standards and calculate the peak area ratios.

$$\text{Peak Area Ratio} = \frac{\text{Moxidectin, Ivermectin or Doramectin Peak Area}}{\text{Abamectin Peak Area}}$$

Construct a linear regression line using the ratios and standard concentrations. The correlation coefficient should be >0.995.

The equation is $y = mx + b$, where

x = Ivermectin, Moxidectin, or Doramectin /Abamectin peak area ratio

y = Ivermectin, Moxidectin, or Doramectin concentration (ppb)

m = slope

b = y-intercept

Incurred tissue ivermectin, moxidectin, or doramectin concentrations should be calculated using this regression line.

Note: Peak heights may be substituted for peak areas if chromatographic peaks display sufficient symmetry.

**United States Department of Agriculture
Food Safety and Inspection Service, Office of Public Health Science**

CLG-AVR.02		Page 12 of 17
Title: Determination of Ivermectin, Doramectin, and Moxidectin by HPLC		
Revision: 02	Replaces: CLG-AVR.01	Effective: 12/22/2006

H. SAFETY INFORMATION AND PRECAUTIONS

1. Required Protective Equipment - Safety glasses, appropriate gloves, lab coat.
2. Hazards

<i>Reagent</i>	<i>Hazard</i>	<i>Recommended Safe Procedures</i>
Acetonitrile, Trifluoroacetic Anhydride, 1-methylimidazole	Flammable and corrosive, may cause skin or respiratory irritation.	Avoid contact or prolonged exposure to vapors. Work in a fume hood. Keep away from flame or heat.
Ivermectin Abamectin Doramectin	Weak teratogen and possible mutagen Severe explosion hazard if in powdered form.	Handle with extreme caution. Handle with extreme caution.
Moxidectin	May cause skin or respiratory irritation. The toxic effects of this material have not been fully evaluated.	Work in a well ventilated area. Store material in a secure, dry, cool well ventilated room.

3. Disposal Procedures

<i>Procedure Step</i>	<i>Hazard</i>	<i>Recommended Safe Procedures</i>
Organic solvents and Avermectin solutions	See above	Collect waste in tightly sealed container and store away from non-compatibles in a cool, well-ventilated, flammable liquid storage area/cabinet for disposal in accordance with local, state, and Federal regulations.

United States Department of Agriculture
Food Safety and Inspection Service, Office of Public Health Science

CLG-AVR.02		Page 13 of 17
Title: Determination of Ivermectin, Doramectin, and Moxidectin by HPLC		
Revision: 02	Replaces: CLG-AVR.01	Effective: 12/22/2006

I. QUALITY ASSURANCE PLAN

1. Performance Standard

<i>Analyte</i>	<i>Analytical Range</i>	<i>Acceptable Recovery (%)</i>	<i>Acceptable Repeatability (CV (%))</i>
Ivermectin	≥ 7.5 ppb	60 -120	< 20
Doramectin	≥ 7.5 ppb	60 -120	< 20
Moxidectin	≥ 7.5 ppb	60 -120	< 20

Acceptability criteria:

- a. Correlation coefficient ≥ 0.995.
- b. Mean recovery for each species in range of 60 -120%.
- c. Repeatability for each species < 20%.
- d. No false positive or false negative results.

2. Critical Control Points and Specifications

<i>Record</i>	<i>Acceptable Control</i>
a. Sample weight	2.5 g ± 0.2 g.
b. Alumina deactivation level	12%
c. 1-methylimidazole/acetonitrile volume	200 µL
d. TFAA/acetonitrile volume	200 µL
e. Condition of SPE column after sample loading (F.1.k.)	Do not allow SPE column to run dry
f. Sample derivatization (F.1.o.)	In dark for a minimum of 15 minutes

3. Readiness To Perform (FSIS Training Plan)

- a. Familiarization
 - i. Phase I: Standards-Duplicate standard curve on each of 3 consecutive days, which will include the following:
 - (a) 0 ppb
 - (b) 7.5 ppb

United States Department of Agriculture
Food Safety and Inspection Service, Office of Public Health Science

CLG-AVR.02		Page 14 of 17
Title: Determination of Ivermectin, Doramectin, and Moxidectin by HPLC		
Revision: 02	Replaces: CLG-AVR.01	Effective: 12/22/2006

- (c) 15 ppb
 - (d) 30 ppb
 - (e) 60 ppb
 - ii. Phase II: Fortified samples - 3 replicates at 0, 7.5, 15, 30 and 60 ppb of both liver and muscle of any applicable species (bovine, ovine, porcine, caprine and equine) over a period of 3 different days.

Note: Phase I and Phase II may be performed concurrently.
 - iii. Phase III: Check samples for analyst accreditation.
 - (a) 8 samples fortified between 0 and 60 ppb. Four of the samples should be liver and 4 muscle. One of the 8 should be blank.
 - (b) Samples submitted by the Quality Assurance Manager (QAM) or supervisor.
 - (c) Letter from QAM is required to commence official analysis.
 - b. Acceptability criteria.

Refer to I.1.
4. Intralaboratory Check Samples
- a. System, minimum contents.
 - i. Frequency: One check sample per week per analyst when samples are analyzed.
 - ii. Records are to be maintained.
 - b. Acceptability criteria.

Refer to I.1.

If unacceptable values are obtained, then:
 - i. Stop all official analyses by that analyst.
 - ii. Take corrective action.
5. Sample Acceptability and Stability
- a. Matrices: Liver, muscle.
 - b. Species: Bovine, ovine, porcine, caprine, and equine.
 - c. Sample size: 16 oz. minimum
 - d. Condition on receipt: Frozen/semi-frozen
 - e. Sample storage:
 - i. Time: 6 months

**United States Department of Agriculture
Food Safety and Inspection Service, Office of Public Health Science**

CLG-AVR.02		Page 15 of 17
Title: Determination of Ivermectin, Doramectin, and Moxidectin by HPLC		
Revision: 02	Replaces: CLG-AVR.01	Effective: 12/22/2006

ii. Condition: Frozen

6. Sample Set

- a. External standard curve at 0, 7.5, 15, 30, and 60 ppb
- b. Recovery
- c. Samples

7. Sensitivity

- a. Minimum proficiency level (MPL): 7.5 ppb.

J. WORKSHEET

The following worksheet is an example.

**United States Department of Agriculture
Food Safety and Inspection Service, Office of Public Health Science**

CLG-AVR.02		Page 16 of 17
Title: Determination of Ivermectin, Doramectin, and Moxidectin by HPLC		
Revision: 02	Replaces: CLG-AVR.01	Effective: 12/22/2006

Worksheet Printed

Avermectin Screen Worksheet

Method

Sample Information	
Analyte Code:	990
Analyst Name:	
Analyst Code:	
Date Started:	
Date Comp'd:	
Check PPB:	
Rec PPB:	

Curve B		Curve A	
Blank	PPB	MOX	ABA
0	7.5		
A Std 1	15		
A Std 2	30		
A Std 3	60		
A Std 4			
Blank	0		
B Std 1	7.5		
B Std 2	15		
B Std 3	30		
B Std 4	60		

Avermectin codes	
Moxidectin:	992
Doramectin:	991
Ivermectin:	923

Sample No	QAType	Curve	MOX	ABA	DOR	IVR	Moxidectin Results		Doramectin Results		Ivermectin Results	
							PPB	%Rec	PPB	%Rec	PPB	%Rec
1												
2												
3												
4												
5												
6												
7												
8												
9												
10												
11												
12												
13												
14												
15												
16												
17												
18												
19												
20												
21												
22												
23												
24												
25												
26												
27												
28												
29												
30												
31												
32												
33												
34												
35												
36												
37												
38												
39												
40												

Regression Output:			
	MOX	DOR	IVR
Y Intercept			
Slope			
R			
No. of Observations	5	5	5
Degrees of Freedom	3	3	3

Regression Output:			
	MOX	DOR	IVR
Y Intercept			
Slope			
R			
No. of Observations	5	5	5
Degrees of Freedom	3	3	3

Analysis Parameters	
500ng/ml Combined Std:	Ivr, Dora, Mox
Combined Std ILN	
500 ng/ml Internal Std:	Avermectin
Internal Std ILN	

With 500 ng / mL conc.:	
Spiking volume of 75uL =	15 ppb
Spiking volume of 150 uL =	30 ppb

Critical Control Points	
Sample weight:	2.5 +/- 0.2g
1-Methylimidazole volume:	200 uL +/- 10 uL
TFAA Volume:	200 uL +/- 10 uL
Final dilution volume:	2.5 ml

LC Parameters	
LC Column(ODS):	15cm X 4.6mm
Mobil Phase:	97:3 MeOH:H2O
Flow Rate:	1.8 ml /min
Iver 1: Injection Volume =	
Iver 2: Injection Volume =	

Blank Used	
Liver ID Number	
Muscle ID Number	

**United States Department of Agriculture
Food Safety and Inspection Service, Office of Public Health Science**

CLG-AVR.02		Page 17 of 17
Title: Determination of Ivermectin, Doramectin, and Moxidectin by HPLC		
Revision: 02	Replaces: CLG-AVR.01	Effective: 12/22/2006

K. APPENDIX

1. References

- a. De Montigny, Pierre, Jung-sook, K. Shim and Pivnichny, J. V., J. Pharm. and Biomedical Anal., Vol. 8 No. 6, pp. 507-511, (1990).
- b. Doherty, Steven J., Fox, Allen and Fink, David W., J.A.O.A.C. Vol. 73, No. 6, pp. 931-934, (1990).
- c. Prabhu, Sunil V., Wherner, Teresa A., Egan, Richard S. and Tway, Patricia C., J. Agric. Food Chem., Vol. 39, pp. 2226-2230, (1991).

L. APPROVALS AND AUTHORITIES

1. Approvals on file.
2. Issuing Authority: Laboratory Quality Assurance Division.