

**United States Department of Agriculture  
Food Safety and Inspection Service, Office of Public Health Science**

SOP No: CLG-TST2.00		Page 1 of 21
Title: Confirmation of Thyreostats by HPLC/MS/MS		
Revision 00	Replaces: NA	Effective: 12/03/04

**Contents**

A.	INTRODUCTION .....	2
B.	EQUIPMENT .....	4
C.	REAGENTS AND SOLUTIONS .....	5
D.	STANDARDS .....	5
E.	SAMPLE PREPARATION .....	6
F.	ANALYTICAL PROCEDURE.....	6
G.	DETECTION AND CONFIRMATION .....	8
H.	SAFETY INFORMATION AND PRECAUTIONS.....	9
I.	QUALITY ASSURANCE PLAN .....	12
J.	WORKSHEET .....	14
K.	APPENDIX .....	15

**United States Department of Agriculture**  
**Food Safety and Inspection Service, Office of Public Health Science**

SOP No: CLG-TST2.00		Page 2 of 21
Title: Confirmation of Thyreostats by HPLC/MS/MS		
Revision 00	Replaces: NA	Effective: 12/03/04

**A. INTRODUCTION**

1. Theory

Thyreostat residues are extracted with acetonitrile from muscle homogenate. The extract is partially cleaned by passing through the silica gel column, and then analyzed by HPLC-MS/MS for confirmation. Confirmation is based on comparison of sample MS/MS spectral data with that of a fortified tissue standard or external standard.







2. Applicability

This method will confirm thyreostats (2-thiouracil, 6-methyl-2-thiouracil, 6-propyl-2-thiouracil, 6-phenyl-2-thiouracil, 2-mercapto-1-methylimidazole, and 2-mercaptobenzimidazole) in porcine, equine and bovine muscle at  $\geq 25$  ppb.

**United States Department of Agriculture**  
**Food Safety and Inspection Service, Office of Public Health Science**

SOP No: CLG-TST2.00		Page 3 of 21
Title: Confirmation of Thyreostats by HPLC/MS/MS		
Revision 00	Replaces: NA	Effective: 12/03/04

3. Structures

THYREOSTATS	
	
2-thiouracil (TU) (MW = 128)	6-methyl-2-thiouracil (MTU) (MW = 142)
	
6-propyl-2-thiouracil (PrTU) (MW = 170)	2-mercapto-1-methylimidazole or Tapazole (TAP) (MW = 114)
	
6-phenyl-2-thiouracil (PhTU) (MW = 204)	2-mercaptobenzimidazole (MBI) (MW = 150)

**United States Department of Agriculture**  
**Food Safety and Inspection Service, Office of Public Health Science**

SOP No: CLG-TST2.00		Page 4 of 21
Title: Confirmation of Thyreostats by HPLC/MS/MS		
Revision 00	Replaces: NA	Effective: 12/03/04

**B. EQUIPMENT**

1. Apparatus

Note: Equivalent apparatus and instrumentation may be substituted for the following:

- a. Robot Coupe® Processor - Robot Coupe U.S.A. Inc.
- b. Centrifuge tubes - 50 mL, polypropylene tube, Falcon Cat. No. 352070, Becton Dickinson Labware.
- c. Pipettors - 5 -100 µL, 100 -1000 µL, Rainin EDP variable volume micropipettes.
- d. Top-loading balance - PM 300, Mettler.
- e. Analytical balance - Leco-250, Leco Corp.
- f. Conical reaction vials - 5mL, Kontes Microflex, #749000-0005.
- g. Microfilterfuge tubes - 45 µm Nylon 66, Cat. No. 7016-22, Rainin.
- h. Evaporator - N-Evap, Organomation Associates.
- i. Nitrogen source - Whatman N2-2010-(75-86).
- j. Volumetric flasks - 1 L, 20, 10 and 1 mL.
- k. Graduated cylinder - 1 L, 500 mL, 10 mL.
- l. HPLC solvent filtering apparatus - with 0.45 µm filter.
- m. Freezer - capable of attaining < - 20 °C.
- n. Pipettes - 0.5 mL, 1 mL and 2 mL.
- o. Centrifuges - Sorvall T6000B and VWR Model V.
- p. Vortex mixer - Genie 2, Fisher Scientific.
- q. Sonicating water bath - Aquasonic, Cat. No. 150T, VWR.
- r. Glass centrifuge tubes - 15 mL, Cat. No. 8084, Pyrex.
- s. Pasteur pipettes - disposable glass, 9 in. long.
- t. HPLC vials and inserts.
- u. SPE column - Silica Gel, Bond Elute, 0.5 g, 10 mL, Cat. No. 1211-3036, Varian.

2. Instrumentation

- a. Micromass Quattro Micro equipped with electrospray LC interface coupled to a Waters 2695 HPLC and autosampler.
- b. LC column - Phenomenex Prodigy 3 µ ODS(3) 100A 4.6 mm x 150 mm.
- c. Guard column - Phenomenex ODS 4 mm x 3 mm.

**United States Department of Agriculture**  
**Food Safety and Inspection Service, Office of Public Health Science**

SOP No: CLG-TST2.00		Page 5 of 21
Title: Confirmation of Thyreostats by HPLC/MS/MS		
Revision 00	Replaces: NA	Effective: 12/03/04

**C. REAGENTS AND SOLUTIONS**

Note: Equivalent reagents and solutions may be substituted.

1. Reagents

- a. Acetonitrile ( $\text{CH}_3\text{CN}$ ) - Cat. No. 015-4, Burdick & Johnson.
- b. Methylene chloride ( $\text{CH}_2\text{Cl}_2$ ) - Cat. No. 9264-03, J. T. Baker.
- c. Methanol ( $\text{MeOH}$ ) - Cat No. 230-4, Burdick & Johnson.
- d. Sodium sulfate - Anhydrous, Cat. No. S421-1, Fisher Scientific.
- e. Formic Acid ( $\text{HCO}_2\text{H}$ ) - Cat. No. 06440, Fluka.
- f. Water - HPLC grade.

2. Solutions

Note: Unless otherwise noted, solutions may be stored at room temperature.

- a. 25%  $\text{MeOH}/\text{CH}_2\text{Cl}_2$  (v/v):  
Mix 1 part  $\text{MeOH}$  with 3 parts  $\text{CH}_2\text{Cl}_2$ .
- b. HPLC mobile phases:  
A = 0.1%  $\text{HCO}_2\text{H}$  - Mix 1 mL  $\text{HCO}_2\text{H}$  with 1L HPLC water.  
B = 0.1%  $\text{HCO}_2\text{H}$  in 1:1  $\text{CH}_3\text{CN}:\text{MeOH}$  - Mix 1 mL  $\text{HCO}_2\text{H}$  with 1 L 1:1  $\text{CH}_3\text{CN}:\text{MeOH}$

**D. STANDARDS**

1. Source

Note: Equivalent sources for the standards can be used.

- 2-Thiouracil - Cat. No. 301507, Aldrich.
- 6-Methyl-2-thiouracil - Cat. No. 69400, Fluka.
- 6-Propyl-2-thiouracil - Cat. No. 82460, Fluka.
- 6-Phenyl-2-thiouracil - Cat. No. P3252, Sigma.
- 2-Mercapto-1-methylimidazole - Cat. No. 301507, Aldrich.
- 2-Mercaptobenzimidazole - Cat. No. M3205, Aldrich.

2. Preparation

- a. Stock Standard Solutions (1 mg/mL):  
Accurately weigh  $20.0 \pm 1$  mg of each of the above standards into separate scintillation vials. Dissolve in 20 mL methanol. Store at  $< -10^\circ\text{C}$ .

**United States Department of Agriculture  
Food Safety and Inspection Service, Office of Public Health Science**

SOP No: CLG-TST2.00		Page 6 of 21
Title: Confirmation of Thyreostats by HPLC/MS/MS		
Revision 00	Replaces: NA	Effective: 12/03/04

- b. Working Standard Solution (3 µg/mL):  
Pipet 30 µL each stock standard to a 10 mL volumetric flask and bring to volume with methanol. Store at room temperature.
- c. LC Standard Solution (25 ng/mL):  
Pipet 10 µL working standard into LC vial and dilute with 400 µL methanol and 800 µL 0.1% formic acid. Prepare as needed and store at room temperature.
- d. Storage and Stability: The stock standard solution is stable for 6 months when stored in a freezer ≤ - 10 °C.

**E. SAMPLE PREPARATION**

After removing excess fat from sample, homogenize with a food processor, transfer into plastic bags and store in a freezer at ≤ - 10 °C. Let the sample thaw prior to analysis.

**F. ANALYTICAL PROCEDURE**

**1. Extraction**

- a. Weigh 5 g homogenized tissue into 50 mL polypropylene centrifuge tube.  
Note: At this time, weigh two 5 g portions of blank muscle tissue into 50 mL polypropylene centrifuge tubes. Use the first tube as a blank and fortify the second tube as a recovery by adding 42 µL of working standard (D.2.b) for a 25 ppb recovery.
- b. Add 10 mL acetonitrile and cap.
- c. Shake vigorously and vortex at least 1 minute until sample is dispersed.
- d. Centrifuge at about 2500 rpm about 5 minutes.
- e. Remove about 5 mL solution and place in 15 mL glass centrifuge tube.
- f. Evaporate on N-Evap to dryness at ≤ 60 °C.
- g. Add 0.5 mL methylene chloride to sample tube, cap and vortex briefly.
- h. Sonicate at least 5 minutes.

**2. SPE Column Cleanup**

Note: Do not let the SPE column to go to dryness during steps 2.b-f below. Also for steps 2.b-f allow liquid level to drain to top of column before adding next volume of liquid.

- a. Add 1 g sodium sulfate to a silica gel SPE column, and position over waste container.
- b. Wash column with 2 mL methylene chloride.
- c. Add sample extract (1.h) to column.

**United States Department of Agriculture**  
**Food Safety and Inspection Service, Office of Public Health Science**

SOP No: CLG-TST2.00		Page 7 of 21
Title: Confirmation of Thyreostats by HPLC/MS/MS		
Revision 00	Replaces: NA	Effective: 12/03/04

- d. Add 0.5 mL methylene chloride to sample tube, vortex and sonicate briefly, and add to column.
- e. Add 2 mL 25% methanol/methylene chloride to sample tube and sonicate at least 1 minute.
- f. Wash column with 2 mL methylene chloride.
- g. Remove waste container and position a 5 mL reaction vial under column.
- h. Elute column with 25% methanol/methylene chloride from (e). Allow eluate to fully drain into tube.
- i. Evaporate to dryness on N-Evap at  $\leq 60^{\circ}\text{C}$ .
- j. Add 200  $\mu\text{L}$  methanol and 400  $\mu\text{L}$  0.1% formic acid.
- k. Cap and vortex briefly.
- l. Add to 0.45  $\mu\text{m}$  filterfuge tubes.
- m. Centrifuge at about 8000 rpm until sufficient volume of filtrate has been collected for HPLC analysis.
- n. Place sample in LC injection vial.

2. Instrument Operating Parameters - LC System

Note: The instrument parameters listed here are examples of one set of suggested optimization parameters. Others may yield equivalent results. The analyst should optimize parameters for the instrument used.

- a. Set initial composition of mobile phase A to 93% and B to 7% at a flow rate of 0.5 mL/min. Allow system to equilibrate.
- b. Injection volume: 20  $\mu\text{L}$
- c. Elution gradient:

Time (min)	Flow Rate (mL/min)	Mobile phase A (%)	Mobile phase B (%)
0	0.5	93	7
6	0.5	93	7
20	0.5	20	80
23	0.5	20	80
25	0.5	93	7
28	0.5	93	7

**United States Department of Agriculture  
Food Safety and Inspection Service, Office of Public Health Science**

SOP No: CLG-TST2.00		Page 8 of 21
Title: Confirmation of Thyreostats by HPLC/MS/MS		
Revision 00	Replaces: NA	Effective: 12/03/04

3. Mass Spectrometer Setup:

Program the mass spectrometer to collect the product ions.

4. Instrumental Settings

Note: Table contains recommended values. Instrumental settings may be adjusted, if necessary, to optimize performance.

Typical LC/MS system setting:

Polarity ES+

Source Temperature 120 °C

Note: See Section K. 1. for additional settings.

5. Injection Sequence

- a. Inject external standard mixture and recovery. Verify that all monitored product ions are present in the external standard.
- b. Inject the recovery and blank. Verify the absence of analyte carry over in the blank. If significant carry over is detected, inject solvent/ blank until reduced to acceptable level.
- c. Inject sample extract(s). Include additional washes and standards in the run as often as necessary to ensure proper identification of sample analytes.
- d. Reinject standard or recovery at the end of the run to verify instrument response.

6. Sample Chromatograms

See Section K.2, Sample Chromatograms

**G. DETECTION AND CONFIRMATION**

1. For each injection:

- a. Plot ion chromatograms for each product ion monitored.
- b. Determine retention times and abundances for all product ions.
- c. Calculate the ratios of product ions specified below for confirming testing.



**United States Department of Agriculture**  
**Food Safety and Inspection Service, Office of Public Health Science**

SOP No: CLG-TST2.00		Page 9 of 21
Title: Confirmation of Thyreostats by HPLC/MS/MS		
Revision 00	Replaces: NA	Effective: 12/03/04

Analyte	Product Ions	Parent mass	Retention Time examples
2-Thiouracil (TU)	84 & 112	129	5.32
2-Mercapto-1-methylimidazole (TAP)	88 & 56	115	6.93
6-Methyl-2-thiouracil (MTU)	126 & 84	143	7.61
6-Propyl-2-thiouracil (PrTU)	154 & 112	171	16.42
2-mercaptobenzimidazole (MBI)	118 & 93	151	17.09
6-Phenyl-2-thiouracil (PhTU)	188 & 146	205	18.63

2. Confirmation Criteria

Confirmation of thyreostat residues in a sample extract requires that the following criteria be met:

- a. Retention time of the product ion peaks in the sample chromatograms must match that found in the external standard or recovery within  $\pm 4\%$ .
- b. At least 2 product ion peaks characteristic of the analyte are present with a signal to noise ratio of greater than 3.
- c. If two product ions are monitored, the presence of one sample ion ratio match that calculated for the external or recovery within a  $\pm 10\%$  arithmetic difference.
- d. The blank does not have confirmable analyte(s).

**H. SAFETY INFORMATION AND PRECAUTIONS**

1. Required Protective Equipment - Protective clothing, eyewear, gloves, and a hood where applicable.

2. Hazards

<i>Procedure Step</i>	<i>Hazard</i>	<i>Recommended Safe Procedures</i>
TU	May be harmful by inhalation, ingestion/skin absorption. Target organs Liver and Thyroid. May cause irritation to skin, mucus membranes and eyes.	Use a hood and wear protective clothing and gloves when handling standards.

**United States Department of Agriculture  
Food Safety and Inspection Service, Office of Public Health Science**

SOP No: CLG-TST2.00		Page 10 of 21
Title: Confirmation of Thyreostats by HPLC/MS/MS		
Revision 00	Replaces: NA	Effective: 12/03/04

MBI	Teratogen. May be harmful by inhalation, in contact with skin and if swallowed. Possible risk of impaired fertility and possible fetus risk.	See above
-----	--	-----------

TAP	Possible Teratogen. May be harmful by inhalation, in contact with skin and if swallowed. Possible risk of impaired fertility and possible fetus risk.	See above
-----	---	-----------

PhTU	May be harmful by inhalation, ingestion and skin absorption. Irritation to mucus membranes and upper respiratory tract.	See above
------	---	-----------

PrTU	Possible Carcinogen. May be harmful by inhalation, ingestion and skin absorption. May cause skin irritation.  Irritation to mucus membranes and upper respiratory tract.	See above
------	--	-----------

MTU	Possible Carcinogen. May be harmful by inhalation, ingestion, and skin absorption. May cause skin irritation. Irritation to mucus membranes and upper respiratory tract.	See above
-----	--	-----------

Acetonitrile/methanol	Flammable, poisonous liquid	Wear protective clothing and gloves when handling acetonitrile and methanol.
-----------------------	-----------------------------	--

**United States Department of Agriculture  
Food Safety and Inspection Service, Office of Public Health Science**

SOP No: CLG-TST2.00		Page 11 of 21
Title: Confirmation of Thyreostats by HPLC/MS/MS		
Revision 00	Replaces: NA	Effective: 12/03/04

Formic acid	Dangerously caustic to skin. Chronic absorption has been reported to cause albuminuria and hematuria.	Wear protective clothing and gloves when handling formic acid.
-------------	---	--

---

3. Disposal Procedures

<i>Procedure Step</i>	<i>Hazard</i>	<i>Recommended Safe Procedures</i>
Thyreostat standards as stated in above table.	See above.	Collect waste in a sealed container and store in a cool, well ventilated, flammable liquid storage area/cabinet for disposal in accordance with local, state, and Federal regulations.
Acetonitrile/methanol	See above.	See above
Formic acid	See above.	See above

**United States Department of Agriculture**  
**Food Safety and Inspection Service, Office of Public Health Science**

SOP No: CLG-TST2.00		Page 12 of 21
Title: Confirmation of Thyreostats by HPLC/MS/MS		
Revision 00	Replaces: NA	Effective: 12/03/04

**I. QUALITY ASSURANCE PLAN**

1. Performance Standard

- a. No false positives from blank tissues.
- b. No false negatives at  $\geq 25$  ppb level fortification.

2. Critical Control Points and Specifications

*Record*

*Acceptable Control*

- a. F.2.b-f. Silica gel column wash      Column must not go to dryness

3. Readiness To Perform (FSIS Training Plan)

a. Familiarization

- i. Phase I: Standards. Inject a mixed standard solution containing all six thyreostats at concentration equivalent to 25 ppb in sample extracts. Repeat analysis on three different days.
- ii. Phase II: Analyst fortified samples. Analyze one blank beef muscle tissue and beef muscle tissue fortified with 25 ppb mixed standards. Repeat the analyses two more days using blank pork muscle for day 2 and blank horse muscle for day 3.

NOTE: Phase I and Phase II may be performed concurrently.

- iii. Phase III: Check samples for analyst accreditation.
  - (a) 6 samples fortified at 25 ppb level of each analyte. Any combination of species may be used, and set must include 1 blank.
  - (b) Report analytical findings to Supervisor and QAM.
  - (c) Letter from QAM is required to commence official analysis.

b. Acceptability criteria.

Refer to I. 1.

4. Intralaboratory Check Samples

a. System, minimum contents.

- i. Frequency: One per week per analyst when samples are analyzed.
- ii. Records are to be maintained.

b. Acceptability criteria.

**United States Department of Agriculture  
Food Safety and Inspection Service, Office of Public Health Science**

SOP No: CLG-TST2.00		Page 13 of 21
Title: Confirmation of Thyreostats by HPLC/MS/MS		
Revision 00	Replaces: NA	Effective: 12/03/04

Refer to I. 1.

If unacceptable values are obtained, then:

- i. Stop all official analyses by that analyst.
- ii. Take corrective action.

5. Sample Acceptability and Stability

- a. Matrices: Bovine, porcine or equine muscles.
- b. Sample receipt, minimum weight: approximately 500 g.
- c. Condition upon receipt: chilled or frozen.
- d. Sample storage:
  - i. Condition: frozen ( $\leq -10$  °C) for blended/homogenized samples.

6. Sample Set

Each set must include the following:

- a. Blank muscle (negative control).
- b. Muscle recovery (positive control).
- c. Samples.

7. Sensitivity

- a. Minimum proficiency level (MPL): 25 ppb

SOP No: CLG-TST2.00		Page 14 of 21
Title: Confirmation of Thyreostats by HPLC/MS/MS		
Revision 00	Replaces: NA	Effective: 12/03/04

## Thyreostat Analysis by LC-MS-MS

Method: TST2.

MeOH/CH<sub>2</sub>Cl<sub>2</sub>:\_\_\_\_\_

[illegible]

**United States Department of Agriculture**  
**Food Safety and Inspection Service, Office of Public Health Science**

SOP No: CLG-TST2.00		Page 15 of 21
Title: Confirmation of Thyreostats by HPLC/MS/MS		
Revision 00	Replaces: NA	Effective: 12/03/04

**K. APPENDIX**

**1. Additional Instrument settings:**

Calibration static	2
Capillary voltage (kV)	4
Cone voltage	35
Extractor voltage	2
RF Lens voltage	0.20
Desolvation Temperature	400 °C
Cone Gas Flow (L/H)	91
Desolvation Gas Flow	512
LM 1 Resolution	15.0
HM 1 Resolution	15.0
Ion Energy	10.4
Entrance	5 - 16
Collision	25 - 21
Exit	1 - 22
LM 2 Resolution	15.0
HM 2 Resolution	15.0
Ion Energy	22.0
Multiplier voltage	22.0
Gas Cell Pirani Pressure (mbar)	5.98e - 3

**United States Department of Agriculture  
Food Safety and Inspection Service, Office of Public Health Science**

SOP No: CLG-TST2.00		Page 16 of 21
Title: Confirmation of Thyreostats by HPLC/MS/MS		
Revision 00	Replaces: NA	Effective: 12/03/04

2. Mass Spectra




Figure 1. Spectra of 25 ppb Thyreostat standards



**United States Department of Agriculture  
Food Safety and Inspection Service, Office of Public Health Science**

SOP No: CLG-TST2.00		Page 17 of 21
Title: Confirmation of Thyreostats by HPLC/MS/MS		
Revision 00	Replaces: NA	Effective: 12/03/04




Figure 2. Spectra of 25 ppb Thyreostats Recovery

**United States Department of Agriculture**  
**Food Safety and Inspection Service, Office of Public Health Science**

SOP No: CLG-TST2.00		Page 18 of 21
Title: Confirmation of Thyreostats by HPLC/MS/MS		
Revision 00	Replaces: NA	Effective: 12/03/04










Figure 3. Spectra of Blank Muscle

**United States Department of Agriculture**  
**Food Safety and Inspection Service, Office of Public Health Science**

SOP No: CLG-TST2.00	Page 19 of 21
Title: Confirmation of Thyreostats by HPLC/MS/MS	
Revision 00	Effective: 12/03/04
Replaces: NA	

3. Possible Fragmentation Patterns of Thyreostats Product Ions

2-thiouracil (TU) (MW = 128)	2-mercapto-1-methylimidazole or tapazole (TAP) (MW = 114)
 <p><math>M+1 = 129^+</math></p>	 <p><math>M+1 = 115^+</math></p>
6-methyl-2-thiouracil (MTU) (MW = 142)	6-propyl-2-thiouracil (PrTU) (MW = 170)
 <p><math>M+1 = 143^+</math></p>	 <p><math>M+1 = 171^+</math></p>
2-mercaptobenzimidazole (MBI) (MW = 150)	6-phenyl-2-thiouracil (PhTU) (MW = 204)
 <p><math>M+1 = 151^+</math></p>	 <p><math>M+1 = 205^+</math></p>

**United States Department of Agriculture**  
**Food Safety and Inspection Service, Office of Public Health Science**

SOP No: CLG-TST2.00		Page 20 of 21
Title: Confirmation of Thyreostats by HPLC/MS/MS		
Revision 00	Replaces: NA	Effective: 12/03/04

3. Reference

Steven J. Lehotay, Alan R. Lightfield, Nichelangelo Anastassiades, and David J. Smith. "Simultaneous Analysis of Beta-Agonists and Thyreostats in Animal Tissues by LC/MS-MS and in-line Fluorescence". 4<sup>th</sup> International Symposium on Hormone and Veterinary Drug Residue Analysis, Antwerp, Belgium, June 4-7, 2002.

**United States Department of Agriculture**  
**Food Safety and Inspection Service, Office of Public Health Science**

SOP No: CLG-TST2.00		Page 21 of 21
Title: Confirmation of Thyreostats by HPLC/MS/MS		
Revision 00	Replaces: NA	Effective: 12/03/04

Approved by:

Date approved:

Eric Flynn

11/29/04

Bill Koscinski

11/29/04

Gina McLeroy

11/29/04

Jess Rajan

11/29/04

Charles Pixley

11/29/04

Phyllis Sparling

11/29/04

Approval signatures on file.